

SECTION SOP-G

EFFECTIVE DATE PAGE

05-01-25 1 OF 1

## **Subject: Table of Contents**

Procedure Number	Subject	Effective Date	Review Date
SOP-G1	Washing of Glassware, Plasticware, and Field Equipment	06-01-25	06-01-25
SOP-G2	Bottle Preparation and Quality Assurance	06-01-25	06-01-25
SOP-G3	Field Analyses and Collection of Samples	Samples are collected by E	•
SOP-G4	Receipt, Handling, and Storage of Samples	06-01-25	06-01-25
SOP-G5	Preparation of Toxicity Samples	06-01-25	06-01-25
SOP-G6	Shipping of Samples	Remo	oved
SOP-G7	Disposal of Samples and Containers	06-01-25	06-01-25
SOP-G8	Reagent Water System Maintenance	06-01-25	06-01-25
SOP-G9	Instrument Maintenance and Repair	06-01-25	06-01-25
SOP-G10	Balance and Weight Calibration	06-01-25	06-01-25
SOP-G11	Mechanical Pipette and Graduated Cylinder Volume Verification	06-01-25	06-01-25
SOP-G12	Thermometer Calibration	06-01-25	06-01-25
SOP-G13	Dessicator Maintenance	06-01-25	06-01-25
SOP-G14	Illumination	06-01-25	06-01-25
SOP-G15	Reagent, Stock Standard, and Chemical Preparation	06-01-25	06-01-25
SOP-G16	Time Calibration	06-01-25	06-01-25



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

4 06-01-25 1 OF 6

SOP-G1

## **Subject: Washing of Glassware and Plasticware**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Julune	06-01-25

## **Document Revision History**

Effective	Revision	Review Type	Evaluators	Revisions
Date	number			
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
12-01-11	1	Internal	Jim Sumner (ETS)	Added quality control requirements for disposable plasticware used in toxicity tests.
09-28-16	2	External	Rick Sherrard (TVA) Don Snodgrass (TVA) Jim Sumner (ETS)	<ul> <li>Updated references to include Standard Methods 22<sup>nd</sup> Edition and NELAP standard.</li> <li>Clarified SOP to indicate that term glassware includes any new or used item in the laboratory, which may include glassware, plasiticware, nitex mesh, field equipment, etc.</li> </ul>
03-01-20	3	Internal	Jim Sumner (ETS)	Updated procedure to include NELAP requirements.
06-01-25	4	Internal	Jim Sumner (ETS)	Updated procedure during document review.



SECTION SOP-G1
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 2 OF 6

### **Subject: Washing of Glassware and Plasticware**

### **Scope and Application**

To clean dirty glassware, plasticware, nitex mesh, and field equipment (The term glassware in this SOP refers to any new or used items in the laboratory, which may include afore mentioned items).

### **Quality Control**

Each batch of items washed must be checked for soap residue using 0.04% bromothymol blue indicator.

### **Equipment and Materials**

Alconox or Sparkleen<sup>™</sup> 1 and 2 detergent powder Nitric acid Sodium bicarbonate Tap water Deionized water Bristle brushes Scouring pads 30-gallon plastic tub Rubber gloves Safety glasses Sink Dishwasher 0.04% Bromothymol blue indicator Storage area for clean glassware Acetone, pesticide grade Rinse bottles Fume hood Detergent Residue Check Log

#### **Procedure**

Note: Gloves and safety glasses must be worn while using nitric acid and acetone.

- A. Preparation of the Nitric Acid Bath (10% nitric acid solution).
  - 1. Add 22.5 L of deionized water to a 30-gal plastic tub.
  - 2. Carefully add 2.5 L of nitric acid to the deionized water and mix.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G1 4 06-01-25 3 OF 6

### **Subject: Washing of Glassware and Plasticware**

#### B. Washing of Glassware and Plasticware.

- New disposable glassware/plasticware do not require thorough cleaning before use. It
  is sufficient to rinse new sample containers with sample or dilution water before use.
  New non-disposable glassware or plasticware must be soaked overnight in 10% nitric
  acid (see procedure below) and rinsed well with deionized water.
  - New disposable sample containers (i.e. cubitainers), medicine cups, and Solo® cups must be checked before use in toxicity tests. Dependent on the use of the consumable, either the survival/reproduction of *Ceriodaphnia* or survival/growth of *Pimephales* are evaluated using the new consumable. Side-by-side comparisons of the organisms with the new and old consumables are used (SOP-AT11 or AT20). Organism survival and reproduction/growth are compared between the old and new consumables. If detrimental effects are noted with the new consumables, they must not be used in toxicity tests.
- 2. All non-disposable glassware/plasticware that comes in contact with samples must be washed after use to remove contaminants in the manner described below.

#### If equipment is washed by hand:

- a. Add a small amount of Alconox or Sparkleen™ 1 detergent (~ 1 Tablespoon) to the sink and fill with hot tap water. Place the dirty glassware in the sink and allow it to soak for at least 15 minutes.
- b. Using scouring pads and bristle brushes, scrub the dirty glassware until it is clean.
- c. Immediately rinse the soapy glassware well with hot tap water.
- d. Test for soap residues using the 0.04% bromothymol blue indicator. Choose a random piece of glassware and place a few drops inside the vessel. If the indicator turns blue, repeat steps a-c. Document the positive/negative result in the Detergent Residue Check Log along with date, analyst initials, and indicator chemical number (Exhibit G1.1).

#### If equipment is machine-washed:

a. Add a small amount of Alconox or Sparkleen<sup>™</sup> 2 detergent to the soap tray in the dish washer.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G1 4 06-01-25 4 OF 6

### **Subject: Washing of Glassware and Plasticware**

- b. Start the dish washer.
- c. After the dish washer's cycle is complete, test for soap residues using the 0.04% bromothymol blue indicator. Choose a random piece of glassware and place a few drops inside the vessel. If the indicator turns blue, repeat steps a-b. Document the positive/negative result in the Detergent Residue Check Log along with date, analyst initials, and indicator chemical number (Exhibit G1.1).

To remove scale and metals:

a. Soak clean glassware in the 10% nitric acid bath. Remove the glassware from the nitric acid bath and rinse well with deionized water.

To remove organic compounds (oily residues):

a. Rinse clean glassware with full-strength, acetone under the fume hood and then triple rinse with deionized water.

Place the clean glassware upside down in the storage area to dry.

Glassware/plasticware for toxicity tests should be thoroughly rinsed with dilution water immediately before use in cultures or toxicity tests.

### C. Disposal of the Nitric Acid Bath.

- 1. Check the nitric acid bath before each use. If the nitric acid becomes dirty or turns yellow, the nitric acid bath must be renewed.
  - a. Sodium bicarbonate is first used to neutralize the acid. Once neutralized, the acid bath may be disposed down the drain with tap water.
  - b. Once completed, see step A.1 for the preparation of the 10% nitric acid bath.

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.



SECTION REVISION NUMBER EFFECTIVE DATE PAGE SOP-G1 4 06-01-25 5 OF 6

### **Subject: Washing of Glassware and Plasticware**

#### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms,  $4^{th}$  ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-014. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.

#### **Exhibits**

Exhibit G1.1: Detergent Residue Check Log





SECTION SOP-G1
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 6 OF 6

## **Subject: Washing of Glassware and Plasticware**

#### **Exhibit G1.1: Detergent Residue Check Log**



Page \_\_\_\_

# Detergent Residue Check Log (0.04% Bromothymol Blue Indicator)

Date	Analyst	Detergent Re	sidue Check
		Positive (BLUE) Rewash Glassware	Negative (YELLOW/GREEN

Note: B = back dishwasher, H = hand washed, No mark = front dishwasher in dish room.

SOP G7-Revision 4-Exhibit G1.1



SECTION SOP-G2
REVISION NUMBER 1
EFFECTIVE DATE 06-01-25
PAGE 1 OF 4

### **Subject: Bottle Preparation and Quality Assurance**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infuse	06-01-25

### **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Type		
06-01-11	0	Internal	Jim Sumner (ETS)	Original document
06-01-25	1	Internal	Jim Sumner (ETS)	Updated procedure during document review.

### **Scope and Application**

To prepare bottles for clients to use for collecting samples.

### **Equipment and Materials**

250-mL plastic bottles
500-mL plastic bottles
1000-mL plastic bottles
Trace Metals Grade Nitric Acid
A.C.S. Certified Sulfuric Acid
Glass Pipettes, pipette bulbs
Permanent red marker
Permanent black marker
Bottle labels



SECTION SOP-G2
REVISION NUMBER 1
EFFECTIVE DATE 06-01-25
PAGE 2 OF 4

### **Subject: Bottle Preparation and Quality Assurance**

#### **Procedure**

#### A. Ordering bottles.

Order unpreserved bottles from an approved vendor.

#### B. Preparation of bottles for chemical analyses

- 1. Preparation of bottles for analyses that require Nitric Acid.
  - a. Clients that are collecting samples for lead and copper must be provided unpreserved 1000-mL bottles. Samples must be received by the laboratory and preserved within 14 days of collection. Upon receipt, each bottle must receive 2.0 mL of nitric acid (to be preserved for analysis). The date that the sample bottle is preserved is recorded on the top of the bottle lid with a permanent marker.
  - b. The amount of nitric acid used for each bottle size is provided in the table below. Use a glass pipette to deliver the required amount of acid to each bottle.

Bottle size (mL)	Amount of Nitric acid (mL)
250	0.5
500	1.0
1000	2.0

c. All bottles that are prepared with Nitric Acid are marked with  $HNO_3$  on the top of the bottle lid with a permanent red marker.



SECTION SOP-G2
REVISION NUMBER 1
EFFECTIVE DATE 06-01-25
PAGE 3 OF 4

### **Subject: Bottle Preparation and Quality Assurance**

- 2. Preparation of bottles for analyses that require Sulfuric Acid.
  - a. The amount of sulfuric acid used for each bottle size is provided in the table below. Use a glass pipette to deliver the required amount of acid to each bottle.

Bottle size (mL)	Amount of Sulfuric acid (mL)
250	0.5
500	1.0
1000	2.0

b. All bottles that are prepared with Sulfuric Acid are marked with a on the top of the bottle lid with a permanent red marker.

#### C. Preparation of sample kits.

- 1. Analysis performed in house.
  - a. Prepare all sample kits in accordance with the client's needs and/or permit requirements. Sample containers are prepared according to test requirements. The volume and preservation requirements for each test performed by the laboratory are provided in QAP Q5: Table Q5.1.
  - b. Place a label on each bottle with the correct preservative listed. This should be written with a permanent black marker.
  - a. Provide chain-of-custody forms with each kit ordered.
  - b. If kits must be shipped, pack coolers to preserve the integrity of the bottles.
- 2. Analysis performed in a sub-contract laboratory.
  - a. Sub-contract laboratories will provide all bottles to Environmental Testing Solutions, Inc. for testing that will be performed by that laboratory.
  - b. Prepare all sample kits in accordance with the client's needs.
  - c. Place a label on each bottle with the correct preservative listed. This should be written with a permanent black marker.



SECTION SOP-G2
REVISION NUMBER 1
EFFECTIVE DATE 06-01-25
PAGE 4 OF 4

### **Subject: Bottle Preparation and Quality Assurance**

- d. Provide chain-of-custody forms with each kit ordered.
- e. If kits must be shipped, pack coolers to preserve the integrity of the bottles.

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 1 OF 15

## **Subject: Receipt, Handling, and Storage of Samples**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Julane	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Type		
12-01-00	0	Internal	Jim Sumer (ETS)	Original document
09-01-09	1	External	William Rogers (TVA)	The sample receipt log, including the assignment of project and sample
		(TVA,	Cynthia Russell (TVA)	numbers, is completed as samples are received in the laboratory.
	Environmental Rick Sherrard (TVA)		Rick Sherrard (TVA)	
	Standard, Inc.) Rock Vitale		Rock Vitale	
	(Environmental		(Environmental	
		Standards, Inc.)		
		Internal	Jim Sumner (ETS)	
08-01-11	2.	External	Lance Ferrell	NC 72-hour hold time exception removed from SOP.
	(NC DENR) (NC DENR)		(NC DENR)	
	Internal Jim Sumner (ETS)		Jim Sumner (ETS)	
01-03-12	3	Internal	Jim Sumner (ETS)	Updated exhibits during document review.
07-01-20	4	Internal	Jim Sumner (ETS)	Procedure for screening toxicity samples for chlorine was moved to SOP-
				C8 and removed from this SOP.
06-01-25	5	Internal	Jim Sumner (ETS)	Updated exhibits during document review.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 2 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

### **Scope and Application**

To provide instructions to ensure the proper receipt, handling, and storage of water samples used for toxicity and analytical tests. This procedure also provides instructions for documenting sample custody.

### **Equipment and Materials**

Mercury-filled or hand-held thermometers
Sample containers
Shipping container (cooler)
Indelible ink pen
Refrigerator
Chain-of-Custody Form
pH Strips
DPD Total Chlorine powder
1-oz Medicine cups
Sample Receipt Log
Sample Condition Form

### Procedure for Receipt, Handling, and Storage of Toxicity Samples

- A. Transfer of Sample Custody for Samples Delivered by Facility Personnel.
  - 1. Upon receipt at the laboratory, a member of the laboratory inspects the shipping container (cooler) to ensure that the container was not tampered with and that the custody seals are intact (if used).
  - 2. The cooler is opened and the sample container (Cubitainer®) is removed and inspected to ensure that the lid remained secure during transport.
  - 3. The chain-of-custody form (Exhibit G4.1) is examined for completeness and correctness of information. At a minimum, the review should ensure that:
    - All pertinent blanks on the form are completed.
    - The sample collection date and times are completed.
    - The sample collector signed and dated the form.
    - The client/facility name, county, permit number, state, and outfall are completed
    - Test requirements (including test dilutions, test species and test type) and sample information are complete.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 3 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

- 4. If the sample container is not labeled, write the name of the sample and the sample date on the Cubitainer® with an indelible ink marker.
- 5. Have the person delivering the sample complete or correct the chain-of-custody form, if necessary, before accepting custody of the sample.
- 6. Measure and record the temperature of the sample (SOP-C1).
- 7. Have the person delivering the sample relinquish custody by printing their name, signing, and recording the date and time in the "Relinquished to ETS by:" line on the form.
- 8. A member of the laboratory then prints their name, signs, and records the date and time in the "Received at ETS by:" line on the form.
- 9. As samples are received, complete the Sample Receipt Log (Exhibit G4.2) by entering the following information into the Excel spreadsheet:
  - Date and time received.
  - Name of laboratory analyst that received the sample.
  - Name of the person that the sample was received by.
  - Sample temperature.
  - Sample name and description.
- 10. As samples are received, assign the sample a unique project and sample number.
  - a. Project numbers are assigned by the next sequential project number for the laboratory. Multiple samples for the same project are given the same project number.
  - b. Sample numbers are assigned with an 8-digit number consisting of the Year, Month, Date, and the sequential number for the sample received by the laboratory on that day (YYMMDD.NN). For example, the second sample received by the laboratory on February 24, 2025 would be assigned the following sample number: 250224.02.
  - c. Enter the project and sample numbers into the sample receipt log and record on chain-of-custody form.
- 11. If requested, make a photocopy of the chain-of-custody form for the person delivering the sample to the laboratory.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 4 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

12. A copy of all chain-of-custody forms is maintained in the project folder in the laboratory files.

#### B. Samples Received by Overnight Carrier.

- 1. Upon receipt at the laboratory, inspect custody seals (if used) to ensure that the shipping container (cooler) and samples were not tampered with. If the custody seals are broken, contact the Laboratory Supervisor to determine appropriate action and document that the custody seals were not intact upon receipt on the chain-of-custody form (Exhibit G4.1).
- 2. After inspection, the cooler is opened and the sample container (Cubitainer®) is removed and inspected to ensure that the lid remained secure during transport.
- 3. Review the chain-of-custody form for completeness and correctness of information. At a minimum, the review should ensure that:
  - All pertinent blanks on the form are completed.
  - The sample collection date and times are completed.
  - The sample collector signed and dated the form.
  - The client/facility name, county, permit number, state, and outfall are completed.
  - Test requirements (including test dilutions, test species and test type) and sample information are complete.
- 4. If the sample container is not labeled, write the name of the sample and the sample date on the Cubitainer® with an indelible ink marker.
- 5. Measure and record the temperature of the sample (SOP-C1).
- 6. Check the appropriate boxes if the custody seals were intact and samples were received in good condition.
- 8. Record the carrier's name and tracking number (if provided) that the sample was delivered by in the "Relinquished to ETS by:" line on the chain-of-custody form.
- 9. A member of the laboratory then prints their name, signs, and records the date and time in the "Received at ETS by:" line on the form.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 5 OF 15

### Subject: Receipt, Handling, and Storage of Samples

- 10. As samples are received, complete the Sample Receipt Log (Exhibit G4.2) by entering the following information into the Excel spreadsheet:
  - Date and time received.
  - Name of laboratory analyst that received the sample.
  - Name of the delivery company that the sample was received by.
  - Sample temperature.
  - Sample name and description.
- 11. As samples are received, assign the sample a unique project and sample number according to A.10.
- 12. If the form is not complete, contact the facility and obtain the required information. Have the facility complete the chain-of-custody form and email a copy to the laboratory. Facility personnel are required by regulatory authorities to correct any errors or omissions they have made to the COC forms. Corrections of facility errors or omissions must not be made by laboratory personnel. Both the incomplete and complete COC forms are included in the test report submitted to the facility.
- 13. A copy of all chain-of-custody forms is maintained in the project folder in the laboratory files.

#### C. Acceptable Temperature upon Receipt

- 1. Samples must be received at the laboratory at 0.0 to 6.0°C. Only grab samples collected relatively short distances from the laboratory may be received above 6.0°C. This exception is only provided if a legitimate effort is made to cool the sample. A legitimate effort is considered if the sample is buried in crushed ice in its shipping container immediately after collection and if the sample arrives at the laboratory within 3 hours of collection. This must be documented on the chain-of-custody form.
- 2. Samples received at the laboratory above 6.0°C apart from the above circumstance will be considered invalid.

#### D. Sample Storage.

1. Samples collected and shipped or delivered to the laboratory must be chilled to 0.0 to 6.0°C and received on ice. Sufficient ice should be placed in the shipping container to ensure that ice will still be present when the sample arrives at the laboratory.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 6 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

- 2. If samples are not to be analyzed upon receipt by the laboratory, they must be stored in a darkened refrigerator at 0.0 to 6.0°C.
- 3. Any unused portions of samples not to be disposed of must be stored in a darkened refrigerator at 0.0 to 6.0°C.

#### E. Sample Holding Times.

1. Each sample must first be used within 36-hours of completion of the sampling period.

#### F. Sample Disposal.

1. All excess sample and test water are disposed of by established laboratory protocols (SOP-G7).

### Procedure for Receipt, Handling, and Storage of Analytical Samples

Note: All samples received by the laboratory are collected by facility personnel and are hand delivered or shipped through overnight delivery.

#### Transfer of Sample Custody for Samples Delivered by Facility Personnel.

- 1. Upon receipt at the laboratory, a member of the laboratory reviews the chain-of-custody form (Exhibit G4.3) for completeness and correctness of information. At a minimum, the review should ensure that:
  - All pertinent blanks on the form are completed.
  - The sample collection date and times are completed.
  - The sample collector signed and dated the form.
  - Test requirements and sample information is complete.
- 2. If the sample container is not labeled, write the name of the sample and the sample date on the bottle with an indelible ink pen or use a bottle label.
- 3. Have the person delivering the sample complete or correct the chain-of-custody form, if necessary, before accepting custody of the sample.
- 4. Measure and record the temperature of the sample (SOP-C1). If applicable, verify that the correct preservations were used (i.e. verify pH requirements using pH strips and verify dichlorination using DPD Total Chlorine powder).



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 7 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

- 5. Have the person delivering the sample relinquish custody by printing their name, signing, and recording the date and time in the "Relinquished by:" line on the form.
- 6. A member of the laboratory then signs their name and records the date and time in the "Accepted by:" line on the form.
- 7. Assign the sample a unique project and sample number and record in the Sample Receipt Log (Exhibit G4.4).
  - a. Project numbers are assigned by the year, date and a sequential number beginning with 500 (YYMMDD.500).
  - b. Sample numbers are assigned with a sequential number.
- 8. Give a copy of the chain-of-custody form (pink copy) to the person delivering the sample to the laboratory.
- 9. A copy of all chain-of-custody forms is maintained in the project folder in the laboratory files.

#### B. Samples Received by Overnight Carrier.

- 1. Upon receipt at the laboratory, open the shipping container. Review the chain-of-custody form (Exhibit G4.3) for completeness and correctness of information. At a minimum, the review should ensure that:
  - All pertinent blanks on the form are completed.
  - The sample collection date and times are completed.
  - The sample collector signed and dated the form.
  - Test requirements and sample information are complete.
- 2. If the sample container is not labeled, write the name of the sample and the sample date on the bottle with an indelible ink pen or bottle label.
- 3. Measure and record the temperature of the sample (SOP-C1). If applicable, verify that the correct preservations were used (i.e. verify pH requirements using pH strips and verify dichlorination using DPD Total Chlorine powder).
- 4. A member of the laboratory then prints their name, signs, and records the date and time in the "Accepted By:" line on the form.
- 5. Assign the sample a unique project and sample number according to A.7.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 8 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

- 6. If the form is not complete, contact the facility and obtain the required information.
- 7. A copy of all chain-of-custody forms is maintained in the project folder in the laboratory files.

#### C. Acceptable Temperature upon Receipt

- Samples must be received at the laboratory at 0.0 to 6.0°C. Only grab samples, collected relatively short distances from the laboratory, may be received above 6.0°C. This exception is only provided if a legitimate effort is made to cool the sample. A legitimate effort is considered if the sample is buried in ice in its shipping container immediately after collection and if the sample arrives at the laboratory and demonstrates a downward trend. The initial field temperature must be documented on the chain-of-custody form.
- 2. If the temperature of the sample is above 6.0°C, and no initial field temperature is provided on the chain-of-custody the client must be notified, and a determination is to be made about recollecting the sample. If another sample cannot be secured a Sample Condition Form (Exhibit G4.5) must be completed and emailed to the NC DEQ regional office.

#### D. Sample Preservation and Dechlorination.

 If any sample that requires preservative or dechlorination (as indicated in QAP-Q5 or Chemistry/Bacteria SOP) is not properly preserved or dechlorinated, the client must be notified, and a determination is made about recollecting the sample. If a sample cannot be secured a Sample Condition Form must be completed and emailed to the NC DEQ regional office (Exhibit G4.5).

#### E. Sample Storage.

- 1. Samples collected and shipped or delivered to the laboratory must be chilled to 0.0 to 6.0°C and received on ice. Sufficient ice should be placed in the shipping container to ensure that ice will still be present when the sample arrives at the laboratory.
- 2. If samples are not to be analyzed upon receipt by the laboratory, they must be stored in a darkened refrigerator at 0.0 to 6.0°C.
- 3. Any unused portions of samples not to be disposed of must be stored in a darkened refrigerator at 0.0 to 6.0°C.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 9 OF 15

### Subject: Receipt, Handling, and Storage of Samples

#### F. Sample Holding Times.

- 1. Please refer to each individual SOP for sample holding times.
- 2. If a sample is out of hold time, the client must be notified, and a determination is made about recollecting the sample. If a sample cannot be secured, then a Sample Condition Form must be completed and emailed to the NC DEQ regional office (Exhibit G4.5).

#### G. Sample Disposal.

1. All excess sample and test water are disposed of by established laboratory protocols (SOP-G7).

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

#### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms, 4<sup>th</sup> ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-014. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 10 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

### **Exhibits**

Exhibit G4.1: Whole Effluent Toxicity Chain-of-Custody Form. Exhibit G4.2: Whole Effluent Toxicity Sample Receipt Log.

Exhibit G4.3: Analytical Chain-of-Custody Form. Exhibit G4.4: Analytical Sample Receipt Log.

Exhibit G4.5: Sample Condition Form.



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 11 OF 15

## **Subject: Receipt, Handling, and Storage of Samples**

### Exhibit G4.1: Whole Effluent Toxicity Chain-of-Custody Form.


Sweeney WTP Purchase order:  Americamysis bahia Effluent dilution: 22.7%	Environmental Testing Solutions, Inc.	nole Effluent Toxicity ain-of-Custody Form	351 Depot Stre Asheville, NC 288 Phone: (828) 350-93 Fax: (828) 350-93
Grab sample: Start date: Time: Volume collected for testing: Composite sample: Start date: Time: Method of transport to laboratory: End date: Time: Method of transport to laboratory: End date: Time: Comments:  End date: Time: Comments:  Chilled during collection (yes/no)? If chilled, specify temperature (°C):  Triple rinse sample container with sample before filling. Completely fill the sample container with no air space. Pack the sample container completely in ice. The sample must be < 6.0°C upon receipt at the laboratory.  Sample custody: (to be completed by sample collector and facility personnel)  BY SIGNING BELOW, I CERTIFY THAT THE PERMIT AND TEST REQUIREMENTS IDENTIFIED ON THIS FORM ARE ACCURATE.  Sample collected by:  Received by:	cility: Sweeney WTP	Purchase order:  Effluent dilution: 22.7%	1 County: New Hanover
Pack the sample container completely in ice. The sample must be < 6.0°C upon receipt at the laboratory.  Imple custody: (to be completed by sample collector and facility personnel)  BY SIGNING BELOW, I CERTIFY THAT THE PERMIT AND TEST REQUIREMENTS IDENTIFIED ON THIS FORM ARE ACCURATE.  Sample collected by:    Print	Grab sample: Start date: Time: Composite sample: Start date: Time: End date: Time: Flow proportional (*): or # samples/hc Chilled during collection (yes/no)?	Sample location:  Volume collected for testing;  Number of containers filled for test  Method of transport to laboratory;  Comments:	
Relinquished by:    Print   Signature   Date and time   Print   Signature   Date and time   Received by:	BY SIGNING BELOW, I CERTIFY THAT THE PERM		FORM ARE ACCURATE.
	Title Agrana		



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 12 OF 15

### **Subject: Receipt, Handling, and Storage of Samples**

#### **Exhibit G4.2: Whole Effluent Toxicity Sample Receipt Log.**





SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 13 OF 15

## **Subject: Receipt, Handling, and Storage of Samples**

### Exhibit G4.3: Analytical Chain-of-Custody Form.

РО Во	omental Testing ox 7565, Ashevill 28) 350-9364, Fa	e, NC 2	8802	1			CH	Cli Ad	ent ddre	Na ess:	ame	8		ST	E			Fax					Upon Recei
Sample 1	dentification	G = Grab, C = Composite	Sample Number ETS Use Only	Matrix: Water, Solid,	Date Collected (MM/DD/YY)	Time Collected (HH:MM)	# Containers	served		),	) DH	NaCH + Dechlor.							1	Field mp. (°c		/	upon Recei
5 5 Sample Condition Temperature (°C) Received on Ice: Sealed Cooler:	Y / N Y / N	inquished	i By	[c	Company		ate			Tire	ne	Acc	epte	d By				Comp	oany		Dat	te	Time
Samples Intact: Additional Comme	Y / N											Prin	nted	Name Name	e of S	amp	re						Date



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 14 OF 15

## **Subject: Receipt, Handling, and Storage of Samples**

### **Exhibit G4.4: Analytical Sample Receipt Log.**

Date received	Client name	Sample identification	Date collected	Time collected	Project number	Sample number	Analyses
27-May-25	Polk Central	Effluent - Grab	27-May-25	0704	250527.500	298145	BOD, 5 day Solids, Total Suspender Fecal Coliform
27-May-25	Adventure Village	Effluent - Composite	27-May-25	0'835	250527.501	298146	BOD, 5 day Solids, Total Suspender Ammonia Nitrogen Fecal Coliform
27-May-25	Saluda WWTP	Influent - Composite	27-May-25	0715	250527.502	298147	BOD, 5 day Solids, Total Suspende
		Effluent - Composite	27-May-25	0709		298148	BOD, 5 day Solids, Total Suspende Ammonia Nitrogen
		Effluent - Grab	27-May-25	0656	The second second	298149	Fecal Coliform
27-May-25	Royal Water Works	Toxaway Falls - Grab	27-May-25	0825	250527.503	298150	BOD, 5 day Solids, Total Suspende Fecal Coliform
27-May-25	Rosman WWTP	Influent - Grab	27-May-25	0830	250527.504	298151	BOD, 5 day Solids, Total Suspende
		Effluent - Composite	27-May-25	0830		298152	BOD, 5 day Solids, Total Suspende
		Effluent - Grab	27-May-25	0830		298153	Fecal Coliform
27-May-25	Sapphire Lakes	Effluent - Grab	27-May-25	0745	250527.505	298154	BOD, 5 day Solids, Total Suspende Fecal Coliform Mixed Liquor Suspende
27-May-25	Columbus WWTP	Influent - Composite	27-May-25	0800	250527.506	298155	BOD, 5 day Solids, Total Suspende
		Effluent - Composite	27-May-25	0805		298156	BOD. 5 day Solids, Total Suspende Ammonia Nitrogen
		Effluent - Grab	27-May-25	0805	The second second	298157	Fecal Coliform
27-May-25	Hot Springs	Influent - Grab	27-May-25	0834	250527.507	298158	BOD, 5 day Solids, Total Suspende
		Effluent - Grab	27-May-25	0832		298159	BOD, 5 day Solids, Total Suspende Ammonia Nitrogen Fecal Coliform
27-May-25	HSHA	Effluent - Grab	27-May-25	0803	250527.508	298160	BOD, 5 day Solids, Total Suspende Fecal Coliform
27-May-25	Baxter	Influent - Composite Effluent - Composite	26-May-25 26-May-25	0802 0810	250527,509	298161 298162	BOD, 5 day BOD, 5 day



SECTION SOP-G4
REVISION NUMBER 5
DATE 06-01-25
PAGE 15 OF 15

### Subject: Receipt, Handling, and Storage of Samples

#### **Exhibit G4.5: Sample Condition Form.**



# SAMPLE CONDITION UPON RECEIPT (SCUR) DEVIATION

15A NCAC 02H .0805 (a) (7) (M): Sample preservation shall be verified and documented. If a laboratory receives a sample subject to G.S. 143-215.1 and 143-215.63 that does not meet sample collection, holding time, or preservation requirements, the laboratory shall document the incident, notify the sample collector or client, and secure another sample that meets the regulatory requirements, if possible. If another viable sample cannot be secured, the original sample may be analyzed but the results reported shall be qualified with the nature of the sample collection, holding time, or preservation infractions and the laboratory shall notify the State Laboratory of the infractions. The notification shall include a statement indicating corrective action taken to prevent future infractions.

The sample(s) identified below deviated from required preservation, hold time, sampling protocol or sample documentation.

Attention:	Mr. Jason Smith (NC DEQ, Laboratory Certification Branch)				
Date:					
Received from:					
Name					
Sample collec	tor				
Phone numbe	r				
Deviation:					
Sample date/time					
Analysis					
Deviation	☐ The sample(s) was not received within holding time.				
	☐ The sample(s) was not received on ice (0 to 6.0°C).				
	☐ The sample(s) was not received within holding time.				
	☐ The sample(s) was not collected and/or preserved correctly (e.g. head space in volatiles, improper container, or preservation/dechlorination)				
	☐ The chain-of-custody did not have all the appropriate information (e.g. collector's name, date collected, time collected, sample identification, number of containers for each analysis).				
	☐ Other:				
Action taken:					
initial delivering	☐ Sample(s) accepted and analyzed per client request.				
	SCUR and COC e-mailed to NC DEQ, Asheville Regional Office – Laboratory Certification Branch.				

NC Wastewater / Groundwater Certificate Number: Wastewater: 600

SOP Q3-Revision 5 – Exhibit Q3.2



SECTION SOP-G5
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 1 OF 6

## **Subject: Preparation of Toxicity Samples**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

### **Document Revision History**

Effective	Revision	Review Type	Evaluators	Revisions
Date	number			
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
07-01-13	1	External	Lance Ferrell	Included clarification that the UV sterilization of samples is only
		(NC DENR)	(NC DENR)	performed for NPDES permits, which includes this testing modification.
				Testing in support of North Carolina NPDES permits are not treated prior to
		Internal	Jim Sumner (ETS)	testing.
03-01-20	2	External (TVA)	Rick Sherrard (TVA)	The use of SSW for NC testing was removed.
				Updated procedure to include NELAP requirements.
		Internal	Jim Sumner (ETS)	Section B.1 corrected to state: Use a 50-mL graduated cylinder to
				measure 10 to 50 mL of sample.
06-01-25	3	Internal	Jim Sumner (ETS)	Updated procedure during document review.



SECTION SOP-G5
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 2 OF 6

### **Subject: Preparation of Toxicity Samples**

### **Scope and Application**

To prepare samples for use in toxicity tests.

### **Equipment and Materials**

60-μm Nitex mesh

Serological, fixed, or adjustable-volume pipettes (Eppendorf®)

Volumetric flasks

Various sized graduated cylinders

Synthetic water

Mercury-filled or hand-held thermometers

Tap water

Temperature controlled incubator

Various sized beakers

Dissolved oxygen (DO) meter equipped with a DO probe

Aeration pumps and aquarium tubing

Serological pipettes

Marine Mix<sup>®</sup> - Forty Fathoms<sup>®</sup> sea salt

Stir bar

Stir plate

Salinity meter

Safety glasses

Gloves

40 watt UV Sterilizer

Various test benchsheets

#### **Procedure**

#### A. Presence of Indigenous Organisms.

1. Filter the water sample through a 60-μm Nitex mesh, if the sample contains indigenous organisms (e.g. copepods or Cladoceran species) that may be confused with the test species. Document on the appropriate toxicity test bench sheet that the sample was filtered through a 60-μm Nitex mesh.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G5 3 06-01-25 3 OF 6

### **Subject: Preparation of Toxicity Samples**

#### B. Preparation of Dilutions.

- If dilutions of the water sample need to be made, use the appropriate volumetric flasks, graduated cylinders and pipettes. The level of precision in the preparation of dilutions is increased when smaller volume graduated cylinders and/or serological pipettes are used to prepare the test concentrations. For example, to measure < 10 mL of sample, use a 10-mL serological pipette. Use a 50-ml graduated cylinder to measure 10 to 50 mL of sample.</li>
- 2. The type of control/dilution water used will depend on the test species (salt or freshwater species).
  - a. For freshwater samples using *Ceriodaphnia dubia, Daphnia* or *Pimephales promelas,* moderately hard synthetic water (MHSW) is used (SOP-AT1). MHSW must have a hardness of 80-100 mg CaCO<sub>3</sub>/L, an alkalinity of 57-64 mg CaCO<sub>3</sub>/L and an initial pH of 6.5-8.5 S.U. (recommended 7.4-7.8 S.U.).
  - b. For marine and estuarine samples using *Americamysis bahia* or *Menidia beryllina*, salt synthetic water (SaltSW) is used (SOP-AT1). SaltSW must have a salinity of 5-32 ppt (typically 25.0  $\pm$  2.0 ppt) and an initial pH of 6.5-8.5 SU.

#### C. Achieving Test Temperature, Equilibration of Dissolved Gases.

- 1. Heat or cool the water samples in a water bath or temperature-controlled incubator, as required, to achieve a final temperature of  $25.0 \pm 1.0^{\circ}$ C ( $26.0 \pm 1.0^{\circ}$ C for *A. bahia* chronic tests)
- 2. Measure the dissolved oxygen (D.O.) concentration (SOP-C2) of each sample after it has been warmed to the required temperature. If the D.O. is > 9.0 mg/L or < 4.0 mg/L, then it may be necessary to aerate the sample. If the sample was warmed too quickly, then the D.O. concentration may be greater than 100% saturation (> 9.0 mg/L).
  - a. Shake the sample to equilibrate the dissolved gases. In most cases, vigorously shaking the sample will lower the D.O. concentration.
  - b. Re-measure the D.O. concentration.
  - c. If the D.O. concentration is still outside the acceptable range, moderately aerate the sample for 1 minute and re-measure the D.O. concentration. Continue aerating for 1-minute intervals and measuring the D.O. concentration until the concentration is within the acceptable range. The sample may be aerated for a total of 15 minutes;



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G5 3 06-01-25 4 OF 6

### **Subject: Preparation of Toxicity Samples**

however, aeration during the preparation of samples should be minimized to prevent the loss of volatile compounds. Document the total aeration time on the appropriate toxicity test bench sheet for the sample.

d. If the D.O. concentration is not within the acceptable range after 15 minutes of aeration, each replicate test vessel must be gently aerated throughout the test period. This activity must be documented on the appropriate toxicity test bench sheet for the sample.

#### D. Saltwater Toxicity Tests, "Salting-Up" the Sample.

- 1. Measure the salinity of the water sample according to SOP-C5. If the salinity is below the required range ( $25.0 \pm 2.0$  ppt unless otherwise specified by the test species and permit requirements), then the sample must be "salted-up" through the addition of sea salts to achieve the required salinity. Consult with the Laboratory Supervisor prior to "salting-up" the sample.
  - a. Pour an aliquot of the water sample into a large beaker with a stir bar.
  - b. Place the beaker on a stir plate and stir the sample.
  - c. Place the tip of the salinity probe in the sample. Measure and record the initial salinity (SOP-C5) on the toxicity test bench sheet for the sample.
  - d. Slowly add Marine  $Mix^{\otimes}$  sea salt to the sample, allow the salinity to stabilize. Continue to add sea salt until the desired salinity is obtained (25.0  $\pm$  2.0 ppt, or as required by the test species and NPDES permit). Record the final salinity on the toxicity test bench sheet for the sample.

#### E. UV Sterilization of Water Samples to Remove Pathogenic Interferences.

*Note:* UV sterilization of samples is only performed for NPDES permits, which include this testing modification.

- 1. To remove pathogenic interferences in minnow chronic toxicity tests, samples may be treated with UV light. A 80-watt UV Sterilizer is used for treating ambient, wastewater, and industrial samples.
  - a. Clean the inside of the UV sterilizer according to procedures outlined in SOP-G1.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G5 3 06-01-25 5 OF 6

### **Subject: Preparation of Toxicity Samples**

- b. Set the UV sterilizer horizontally in the supports on the counter, with the inlet and outlet tubes facing up.
- c. Rinse the inside chamber of the sterilizer with deionized water. Carefully pour approximately 2000 ml deionized water into the outlet tube of the sterilizer. Using a rubber stopper, plug the outlet tube.
- d. Remove the UV sterilizer from the counter, making sure to keep the sterilizer horizontal. Gently rock the sterilizer back and forth to thoroughly rinse the inside chamber.
- e. Remove the rubber stopper and drain the sterilizer.
- f. Beginning with the lowest dilution or control sample, carefully pour the required volume to be treated into the outlet tube. Using a rubber stopper, plug the outlet tube.
- g. Plug the ballast cord into the electrical outlet. Ensure the UV light is illuminated, by looking through the view port.
- h. Set the timer for 2 minutes and turn the timer on.
- Remove the UV sterilizer from the counter, making sure to keep the sterilizer horizontal. Gently rock the sterilizer back and forth to provide the maximum exposure of the sample to the UV light.
- j. Return the UV sterilizer to the supports on the counter, horizontally with the inlet and outlet tubes facing up.
- k. After 1 minute has elapsed, repeat steps i and j.
- I. After 2 minutes have elapsed, unplug the ballast cord from the electrical outlet and turn off the timer.
- m. Remove the stopper from the outlet tube and carefully pour the sample into a clean glass beaker or flask. The UV treated sample must never come in contact with glassware or equipment that has been exposed to non-treated sample.
- n. Repeat steps c through m.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G5 3 06-01-25 6 OF 6

### **Subject: Preparation of Toxicity Samples**

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

#### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms,  $4^{th}$  ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.



SECTION SOP-G7
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 1 OF 2

### **Subject: Disposal of Samples and Containers**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Julane	06-01-25

### **Document Revision History**

Effective	Revision	Review Type	Evaluators	Revisions
Date	number			
05-06-01	0	Internal	Jim Sumner (ETS)	Original document
03-01-20	1	External (TVA)	Rick Sherrard (TVA)	Updated procedure to include NELAP requirements. Corrected grammatical errors.
		Internal	Jim Sumner (ETS)	For TVA sample containers, obliterate the name "TVA" (if present) with a sharpie.
06-01-25	2	Internal	Jim Sumner (ETS)	Updated procedure to include requirements for disposing sample containing tritium and colilert samples.     Removed requirement for TVA sample containers (TVA is not documented on container).

## **Scope and Application**

To properly dispose of samples and containers.

### **Equipment and Materials**

Sample containers (Glass or Plastic) Samples Recycling bin Sharpie

### **Procedure**

#### A. Disposal of Non-hazardous Samples

- 1. Only samples can be disposed of after all analyses have been performed.
- 2. All sample containers must be properly disposed.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G7 2 06-01-25 2 OF 2

### **Subject: Disposal of Samples and Containers**

- a. To purge a sample, turn on the cold water tap and slowly pour the sample down the drain.
- b. After the sample container is empty, rinse it with tap water.
- c. Place the container in a recycling bin.
- d. Recycled materials are taken to the recycling center every week.

#### B. Disposal of Hazardous Samples and Containers.

- 1. All containers that contain hazardous waste must be returned to the client.
- 2. The laboratory must not dispose of any samples that contain asbestos, PCB's or oil contaminated waste.
- 3. Colilert samples are disposed of according to procedures described in SOP-B6.
- 4. Samples containing tritium are disposed of according to procedures described in Policy-P6.

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements. Additional requirements for the disposal of tritiated samples is contained within this policy.



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

3 06-01-25 1 OF 5

SOP-G8

### **Subject: Reagent Water System Maintenance**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	~	06-01-25
Quality Assurance Officer	Jim Sumner	Infunse	06-01-25

### **Document Revision History**

Effective	Revision	Review Type	Evaluators	Revisions
Date	number			
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
06-29-09	1	Internal	Jim Sumner (ETS)	Included a flow diagram for the deionized water system.
03-01-20	2	External (TVA)	Rick Sherrard (TVA)	The Milli-Q water system was removed.
				Increased deionized water system tanks to 3.6 CF.
		Internal	Jim Sumner (ETS)	Updated procedure to include NELAP requirements.
06-01-25	3	Internal	Jim Sumner (ETS)	Updated chlorine requirements.
				Updated reference to 24 <sup>th</sup> Edition of Standard Methods.

### **Scope and Application**

To maintain the deionized water system for cultures and toxicity/analytical tests.

### **Quality Control**

- 1. Total residual chlorine is monitored in deionized water through the testing of blanks (SOP-C8). Blanks are tested with each batch of toxicity samples analyzed for chlorine. Levels of chlorine should be less than  $\frac{1}{2}$  the detection limit (<25  $\mu$ g/L). If detectable results are obtained, corrective action must be taken.
- 2. The activated carbon tank is replaced every 6 months.
- 3. The prefilters and post filters in the deionized water system are replaced every 6 months.
- 4. The working and polishing deionizers are replaced when the resistivity light goes out or every 6 months (whichever comes first).
- 5. All maintenance activities must be recorded in the Reagent Water System Log (Exhibit G8.1).

3

2 OF 5



**SECTION** SOP-G8 **REVISION NUMBER EFFECTIVE DATE** 06-01-25 PAGE

### **Subject: Reagent Water System Maintenance**

### **Equipment and Materials**

Potable water Pressure regulators Course filter Prefilters and housings Prefilter Cartidge, 1 micron Prefilter housing assembly, 20" Activated carbon tank, 3.6 cf Mixed bed deionizers, 3.6 cf Post filter Cartridges, 0.2 micron, 10" Post filter housing assemblies, 10" Interconnecting assembly (includes resistivity lights, pressure regulator, standard FDA grade tubing and fittings, and water meter) Water valve Reagent water system log

#### **Procedure**

#### A. Flow Chart.

- 1. City of Asheville, NC potable water is diverted from the city supply line into the laboratory. The main line passes through a pressure regulator and then through a course filter and prefilter to supply the laboratory with potable water (Figure G8.1).
- 2. Standard FDA grade tubing is diverted from the source potable water line to a pressure regulator.
- 3. From the pressure regulator, water is passed through a 20-inch, 1-micron prefilter and then a 3.6 cubic-foot activated carbon tank.
- 4. From the activated carbon tank, water is passed through a 3.6 cubic-foot working deionizer. After the working deionizer, a resistivity light indicates the quality of megaohm-cm water. The light is illuminated green with high quality water. When the working deionizer is exhausted, the light will turn red. At this time, the working deionizer and activated carbon tank must be replaced.
- 5. From the working deionizer, the water passes through a second 3.6 cubic-foot polishing deionizer.




SECTION SOP-G8
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 3 OF 5

### **Subject: Reagent Water System Maintenance**

- 6. After the polishing deionizer, the water passes through a 10-inch, 0.2-micron post filter.
- 7. From the post filter, deionized water is available.

Figure G8.1: Diagram of the deionized water system.



### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms, 4<sup>th</sup> ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.



SECTION SOP-G8
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 4 OF 5

## **Subject: Reagent Water System Maintenance**

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-014. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.

### **Exhibits**

Exhibit G8.1: Reagent Water System Log.



3



**SECTION** SOP-G8 **REVISION NUMBER EFFECTIVE DATE** 06-01-25 5 OF 5 PAGE

## **Subject: Reagent Water System Maintenance**

### **Exhibit G8.1: Reagent Water System Log.**



Page	

### Reagent Water System Log

**Location:** □ Main Laboratory □ Fish Culture Laboratory

Service Date Analyst	Activity (type of service performed and comments)

SOP G8-Revision 3-Exhibit G8.1



SECTION SOP-G9
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 1 OF 2

### **Subject: Instrument Maintenance and Repair**

### **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

### **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Туре		
09-01-09	0	Internal	Jim Sumner (ETS)	Original document
06-20-12	1	External (TVA)	William Rogers (TVA) Donald Snodgrass (TVA) Rick Sherrard (TVA)	Exhibit G9.1 amended to include the dates of routine maintenance, removed from service and returned to service.
22 24 22		Internal	Jim Sumner (ETS)	
03-01-20	2	Internal	Jim Sumner (ETS)	Updated procedure to include NELAP requirements.
06-01-25	3	Internal	Jim Sumner (ETS)	Updated Instrument Maintenance and Equipment Maintenance
				Logs.
				Updated procedure regarding individual item repairs.

## **Scope and Application**

To provide instruction and documentation on the maintenance and repair of laboratory instruments.

### **Procedure**

### A. Maintenance, Repair, and Nonconformances

- Service contracts or in-house preventive maintenance and/or repair are documented in a central maintenance log, which may include documentation for each instrument (example provided in Exhibit G9.1). Larger equipment (e.g. incubators) repairs may also be recorded on the side of the equipment. Documentation may include annual maintenance and maintenance performed at other times due to malfunction (except for routine cleaning).
- 2. Nonconformances (e.g. exceeding quality control limits) for each instrument are also documented in the maintenance log. This includes what led to the nonconformance, probable cause, and how the non-conformance was resolved. In addition,



SECTION SOP-G9
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 2 OF 2

## **Subject: Instrument Maintenance and Repair**

nonconformances are documented through QAP-Q3, Surveillance and Corrective Action Reports.

- 3. Manuals provided by the instrument manufacturers are maintained on file as reference tools to trouble shoot problems with instrumentation. If problems cannot be resolved internally, service technicians are hired to repair or calibrate the malfunctioning instruments.
- 4. Calibration documentation for each instrument is maintained in various logbooks for the tests performed on that instrument or in calibration logbooks maintained for that instrument. Calibration information is not documented in the maintenance log.

### References

**Instrument Manuals** 

### **Exhibits**

Exhibit G9.1: Example Maintenance Log (attached)



SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 1 OF 8

## **Subject: Balance and Weight Calibration**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Туре		
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
09-01-09	1	External (TVA, Environmental Standard, Inc.)	William Rogers (TVA) Cynthia Russell (TVA) Rick Sherrard (TVA) Rock Vitale (Environmental Standards, Inc.) Jim Sumner (ETS)	Balance calibration range changed to bracket all potential measurements performed in the laboratory. Removed warning limits from acceptance criteria. Corrective action included if certified weights exceed acceptance criteria. Changed balance logs and exhibits to reflect these changes and to provide traceability to the last balance and NIST calibrations performed by external vendors.  NIST certified weights recalibrated, which changed the true values.
01-03-12	2	Internal	Jim Sumner (ETS)	Amended associated exhibit and values throughout document.  • Ammended wall chart to provide the source (SOP G10 – Table G10.1, revision 09-01-09).  • Included calibration procedure for the Cahn 28 Automatic Electrobalance.
				Updated exhibits during document review.
07-01-20	3	Internal	Teresa Sleeper (ETS) Jim Sumner (ETS)	Replaced calibration instructions for Fisher Scientific ACCU-224 with instructions for Mettler-Toledo ME204.      Updated procedure to include NELAP requirements.      Additional guidance included in SOP.      Updated log to reflect new control limits for certified weights.
06-01-25	4	Internal	Jim Sumner	Updated exhibits during document review.



SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 2 OF 8

### **Subject: Balance and Weight Calibration**

### **Scope and Application Purpose**

To maintain balance and weights for toxicity and analytical tests.

### **Quality Control**

- 1. Balances must be serviced yearly by an approved company (e.g. Laboratory Instruments, Inc.). Documentation is maintained in the laboratory's QC files.
- 2. Balances must be verified and/or calibrated before each use.
  - Mettler-Toledo ME204:
    - Internal calibration is performed at least weekly.
    - o Balance level verified and corrected as needed before each use.
    - Balance tared before each use.
    - o Balance verified with at least two Class "S" Weights each day prior to use. The weights must bracket the items to be weighed.
    - o Balance verified with seven Class "S" Weights weekly.
  - Cahn 28 Automatic Electrobalance:
    - o Balance tared and calibrated before each use.
    - Balance verified with two Class "S" Weights before and after use. The weights must bracket the items weighed.
- 2. Class "S" Weights must be verified every 5 years against NIST certified weights. Precision Weighing, Inc. in Cary, NC is used to perform this verification. Documentation is maintained in the laboratory's QC files.
- 3. All certified weight measurements <u>must</u> be within the control limits before any materials or samples can be analyzed. If the measurements are out of range, perform the balance internal calibration and document in the balance log. Re-analyze the certified weights. If the measurements still exceed the control limits, notify the Laboratory Supervisor to determine if the balance requires cleaning or adjustment. If the balance cannot be corrected to meet the control limits through cleaning or adjustment, the balance must be serviced by an approved instrument calibration company (e.g. Laboratory Instruments, Inc.). If final measurements using the Cahn Electrobalance exceed the control limits, all weight measurements performed that day must be re-analyzed. Refer to Table G10.1 for the lower and upper control limits for weights used to calibrate balances.



SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 3 OF 8

## **Subject: Balance and Weight Calibration**

**Table G10.1:** Lower and upper control limits for Class "S" Weights used to calibrate balances.

True value (g)	Control Limit (g) (0.01005 to 0.1005 g weights ± 2.0%) (1.0000 to 149.9998 g weights ± 0.10%)						
	Lower	Lower Upper					
0.01005	0.0985 (9.85 mg) 0.01025 (10.25 mg)						
0.05004	0.04904 (49.04 mg) 0.05104 (51.04 mg)						
0.1005	0.0985	0.1025					
1.0000	0.9990	1.0010					
10.0001	9.9901	10.0101					
49.9999	49.9499 50.0499						
99.9999	99.8999	100.0999					
149.9998	149.8498	150.1498					

## **Equipment and Materials**

Balance (Mettler-Toledo ME204, Cahn 28 Automatic Electrobalance) Class "S" Weights

0.01005, 0.05004, 0.1005, 0.20001, 1.0000, 10.0001, 49.9999, 99.9999 and 149.9998 g

Anti-static brush

Forceps Spatula

Weigh boat or medicine cup

Balance Log



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

SOP-G10 4 06-01-25 4 OF 8

**Subject: Balance and Weight Calibration** 

### **Procedure**

Note: Always use forceps to place weights on the balance.

#### A. Calibration of the Mettler-Toledo ME204.

- 1. The balance must be verified before each use. Record calibration information in the Balance Log (Exhibit G10.1).
- 2. Turn the balance on by pressing the  $\rightarrow 0/T \leftarrow$  button.
- 3. Level the balance using the leveling feet so that the air bubble is centered within the circle of the level indicator and document in the Balance Log.
- 4. Wait until the balance reads 0.0000 g before starting. Use the →0/T← key to tare the reading.
- 5. For weighing chemicals and materials, use two certified weights that bracket the weight of chemicals or materials to be measured. Record the measurement of both weights in the Balance Log.
- 6. Weekly perform the internal calibration of the balance and record in the Balance Log. Press and hold the Cal button until ADJUST appears. Select ADJ.INT by pressing Menu and allow the balance to perform the internal calibration. Once the calibration is complete, the balance will read ADJ.DONE and 0.0000 g will appear when the balance is ready.
- 7. Weekly (when the internal calibration is performed), verify the balance with all the Class "S" Weights identified on the log. Record the measurements in the Balance Log.

#### B. Weighing Chemicals or Materials using the Mettler Toledo ME204 Balance

- Before weighing any materials or at a minimum daily, verify the balance according to B.4

   B.9.
- 2. Place a weigh boat on the balance pan and press  $\rightarrow 0/T \leftarrow$  to tare to 0.0000 g.
- 3. Using a spatula, carefully spoon the chemical or place the material onto the center in a weigh boat. Make sure to close the glass doors to ensure a correct reading. Allow the reading to stabilize and record the measurement.

SOP-G10

06-01-25

5 OF 8



SECTION
REVISION NUMBER
EFFECTIVE DATE
PAGE

## **Subject: Balance and Weight Calibration**

- 4. Once completed, remove the weigh boat and sweep out the interior with an anti-static brush.
- Place the balance in standby by pressing and holding the  $\rightarrow 0/T \leftarrow$  button until **STANDBY** appears. The balance will display **MT.GREEN** while in standby mode.

#### C. Calibration of the Cahn 28 Automatic Electrobalance.

- 1. The balance must be verified before each use. Record calibration information in the Balance Log (Exhibit G10.1).
- 2. Turn the balance on by pressing the **POWER** button (located in the back of the balance). The balance should warm for a minimum of 1-hour before use.
- 3. Set the balance to the A Range (200 mg 1  $\mu$ g).
- 4. Press the **TARE** button on the balance (the balance must read 0.00 mg before starting).
- 5. Place the 200.01 mg weight on balance Pan A and wait until ready. Record the answer in the Balance Log. Press the **CAL** button to correct the reading and record the corrected value in the Balance Log.
- 6. Remove the weight from Pan A and allow the balance to stabilize to zero. If the reading has drifted, recalibrate the balance.
- 7. Place the 10.05 mg weight on balance Pan A and wait until ready. Record the answer in the Balance Log.
- 8. Place the 50.04 mg weight on balance Pan A and wait until ready. Record the answer in the Balance Log.

### D. Weighing Materials (weigh pans used in toxicity tests).

- 1. Before weighing any materials, calibrate the balance according to A.4 A.9.
- Using the forceps, carefully place the weigh pan onto the center of Pan A. Make sure to
  close the glass door to ensure a correct reading. Allow the reading to stabilize and
  record the measurement. Record the date that the measurements are being made and
  your initials on the bench sheet.
- 3. Once completed, remove the pan.



SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 6 OF 8

## **Subject: Balance and Weight Calibration**

- 4. Once all measurements are completed, measure the 10.05 mg and 50.04 mg weights. Allow the readings to stabilize and record the answers in the Balance Log.
- 5. Turn the balance off by pressing the **POWER** button (located in the back of the balance).

## References

Instrument Manual

### **Exhibits**

Exhibit G10.1: Balance Log.





SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 7 OF 8

## **Subject: Balance and Weight Calibration**

### Exhibit G10.1: Balance Log (Mettler-Toledo ME204).



Page \_\_\_\_

### Mettler-Toledo ME204 Balance Log

Instrument	Serial number	Calibration date	Calibration company
Mettler-Toledo ME204 Balance	C006982938	01-16-25	Laboratory Instrument Services
Certified Weights	20410	03-03-23	Precision Weighing

	Internal	Level	Certified weight (g) and control limits*					Analy			
Date	Cal.	(v)	0.0102	0.1005	1.0000	10.0001	49.9999	99,9999	149.9998	initials	
	(4)		0.0100-0.0104	0.0985-0.1025	0.9990-1.0010	9.9901-10.0101	49.9499-50.0499	99.8999-100.0999	149.8498-150.1498	-	
			-								
			-								
			1		1 - 1		100				

\*Certified weights must be within established control limits before measurements are made.

If certified weights exceed control limits, perform balance internal calibration and re-measure certified weights.

SOP G10-Revision 4-Exhibit G10.1



SECTION SOP-G10
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 8 OF 8

## **Subject: Balance and Weight Calibration**

### Exhibit G10.1: Balance Log (Cahn 28 Automatic Electrobalance).



Page \_\_\_\_

### Cahn 28 Automatic Electrobalance Log

Instrument	Serial number	Calibration date	Calibration company
Cahn 28 Automatic Electrobalance	41520	01-16-25	Laboratory Instrument Services
Certified Weights (10.05 and 50.04 mg)	20410	03-03-23	Precision Weighing
Certified Weight (200.01 mg)	1000055140	03-03-23	Precision Weighing

Date	Data	Zero	Zero Tare		on weight 01 mg)	(m	ation weights  ng)  ments are taken)	weigh	rification ts (mg) nents are taken)	Analys
	(v)	Initial	Final (calibrated to 200 mg)	10.05 Centrol Limit* (9.85-10.25)	50.04 Control Limit* (49.04-51.04)	10.05 Control Limit* (9.85-10.25)	50.04 Control Limit* (49.04-51.04)	initials		
	-									
	-									

<sup>\*</sup>Initial and final verification weights must be within established control limits. If verification weights exceed control limits, re-zero and re-calibrate balance and re-analyze samples.

SOP G10-Revision 3- Exhibit G10.1



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 1 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	2	06-01-25
Quality Assurance Officer	Jim Sumner	Jan/unse-	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators Revisions				
Date	number	Type					
12-01-00	0	Internal	Jim Sumner (ETS)	n Sumner (ETS) Original document			
01-03-12	1	Internal	Jim Sumner (ETS)	Updated exhibits during document review. Added leak test to procedure.			
06-20-12	2	External	William Rogers (TVA)	Graduated cylinder volume verification procedure added.			
		(TVA)	Donald Snodgrass (TVA)				
			Rick Sherrard (TVA)				
		Internal	Jim Sumner (ETS)				
10-31-12	3	External	William Rogers (TVA)	Clarified graduated cylinder volume verifications for only polypropylene			
		(TVA)	Donald Snodgrass (TVA)	cylinders.			
			Rick Sherrard (TVA)	Updated exhibit to provide acceptable ranges for polypropylene cylinders			
				and included a table with these limits.			
		Internal	Jim Sumner (ETS)				
06-01-25	4	Internal	Blake Lancaster (ETS)	Updated Z factor table to be based on only current temperature and air			
				pressure			
				Updated pipette and graduated cylinder verification log sheets to include			
				serial numbers and verification start and end times.			
				Updated recommended limits.			

## **Scope and Application**

To maintain the accuracy of mechanical pipettes and graduated cylinders used for toxicity and analytical testing.



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 2 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

## **Quality Control**

- 1. All mechanical pipettes must be calibrated once every 6 months. In addition, mechanical pipettes must pass a leak check. 100 mL and 250 mL polypropylene graduated cylinders must be verified annually.
- 2. All pipettes must be within the recommended limits established by the manufacturer. If the readings are out of range, notify the Laboratory Supervisor and take the pipette out of use. If the pipette cannot be corrected to meet the recommended limits through cleaning or adjustment, the pipette must be sent to an approved instrument calibration company (e.g. Precision Weighing) to be recalibrated. Refer to Table G11.1 for the recommended limits for pipette volumes used to calibrate mechanical pipettes.

Table G11.1: Recommended limits for mechanical pipette volumes

### Repeater Pipettes:

True value (ml)	Recommended Limits (g)		Recommended CV (%)
	Lower	Upper	
0.0500	0.0491	0.0509	< 1.5%
0.1000	0.0982	0.1018	< 1.5%

#### Adjustable Finnpipettes:

True value (ml)	Recommended Limits (g)		Recommended CV (%)
	Lower	Upper	
0.0100	0.0092	0.0108	< 3.0%
0.0500	0.0492	0.0508	< 1.5%
0.1000	0.0992	0.1008	< 1.5%
0.2500	0.242	0.258	< 1.5%
0.5000	0.492	0.508	< 1.5%
1.0000	0.992	1.008	< 1.5%



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 3 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

3. Polypropylene graduated cylinders must be within ± 5% of the true value. If the readings are out of range, notify the Laboratory Supervisor and take the cylinder out of use. Refer to Table G11.2 for the recommended limits for cylinder volumes used to verify polypropylene graduated cylinders.

**Table G11.2:** Recommended limits for polypropylene graduated cylinder volumes

True value (mL)	Recommended Limits (g)			
	Lower	Upper		
100	95 105			

## **Equipment and Materials**

Balance (Mettler-Toledo ME204, or equivalent)
Class "S" Weights
100-mLbeaker
Deionized water
Thermometer
Adjustable-Fixed Mechanical Pipettes (Eppendorf® and Fisherbrand HandyStep®)
Mechanical pipette tips
100 mL and 250 mL Graduated Cylinders
Pipette Volume Verification Logsheet
Graduated Cylinder Volume Verification Logsheet

### **Procedure**

- A. Mechanical Pipette Verification Procedure.
  - 1. Calibrate the balance before pipette verification (SOP-G10).
  - 2. Complete the following on the Pipette Volume Verification Logsheet (Exhibit G11.1):
    - Analyst's initials
    - Date checked
    - Pipette serial number
    - Pipette manufacturer



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 4 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

- Water and ambient temperature
- Volumes checked
- Time started
- 3. Fill a 100-ml beaker with deionized water. Allow the water to adjust to room temperature. Record the temperature on the Logsheet.
- 4. Perform a leak test by placing a tip on the pipette and filling with deionized water to the maximum nominal volume and placing the pipette on a vibration-free stand. Observe the meniscus in the tip for 1 minute. After 1 minute, there should be no visible droplet formation at the tip. If there are no signs of droplet formation, proceed with the calibration check. If droplet formation is observed, successful maintenance and corrective action must be taken or the pipette must not be used.
- 5. Place a weigh boat on the balance pan and **TARE** to 0.0000 mg
- 6. Draw up the required volume of deionized water.
- 7. Release the water into the weigh boat and wait until the balance is stable. Record the reading on the Pipette Volume Verification Logsheet (Exhibit G11.1). Repeat this process 10 times. Record verification end time when completed.
- 8. Continue to verify and record the different pipette volumes.
- 9. Enter the results into the Fixed-Volume Pipette Calibration excel spreadsheet (Exhibit G11.2). Results are corrected for temperature and air pressure, where a Z-factor is applied (Exhibit G11.3).

### B. Polypropylene Verification Procedure.

- 1. Calibrate the balance before graduated cylinder verification (SOP-G10).
- 1. Complete the following on the Polypropylene Graduated Cylinder Volume Verification Logsheet (Exhibit G11.4):
  - Date
  - Analyst
  - Water and ambient temperature
  - Volume verified
  - Cylinder type (100 mL, 250 mL, etc.)



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 5 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

- Cylinder identification number
- 2. Fill an Erlenmeyer flask with deionized water. Allow the water to adjust to room temperature. Record the temperature on the logsheet.
- 3. Place the cylinder on the balance pan and **TARE** to 0.0000 mg
- 4. Fill the cylinder with deionized water to the appropriate volume to verify. Wait until the balance is stable. Record the reading on the Graduated Cylinder Volume Verification Logsheet (Exhibit G11.4).
- 5. Continue to verify and record for each cylinder to be checked.
- 6. Enter the results into the Fixed-Volume Pipette Calibration excel spreadsheet (Exhibit G11.5). Results are corrected for temperature and air pressure, where a Z-factor is applied (Exhibit G11.3).

### References

Instrument Manual

#### **Exhibits**

Exhibit G11.1:	Example Pipette Volume	Verification Logsheet.

Exhibit G11.2: Fixed-Volume Pipette Calibration excel spreadsheet

Exhibit G11.3: Z-factor Correction Factor Table

Exhibit G11.4: Polypropylene Graduated Cylinder Volume Verification Logsheet.

Exhibit G11.5: Polypropylene Graduated Cylinder Volume Verification excel spreadsheet



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 6 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

### **Exhibit G11.1: Example Pipette Volume Verification Logsheet**



Page 1 of 1

### **Pipette Volume Verification**

1		
	End time	
Fisherbrand Handy Step	Serial number	17E59354
	Digital thermometer SN	
Mettler-Toledo ME204	Serial number	C006982938
		Fisherbrand Handy Step Serial number Digital thermometer SN

9		
10		
1-minute Leak Check	Pass 🗆	Fail 🗆
Start time	End time	

SOP G11 - Revision 4 - Exhibit G11.1



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 7 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

### Exhibit G11.2: Fixed-Volume Pipette Calibration excel spreadsheet

				•	Calibration			
			кете	rence: ISO 8655-6:2	2002 (E)			I
Analyst:		JS			Test Conditions			Norma
Date of Calibr	ation:	12-31-2	4		Barometric pres	sure. kPa or	mbar	101.2
Calibration d		06-25	·		Relative Humidi			87
Pipet manufa		Finnpip	ette		Ambient temp., o			23.0
Serial/ID num		U98294			Test liquid temp			23.0
Nominal Volu		0.01			Z correction fac		1.0035	•
Location of us		Chemis	ry lab					
eak Test					Balance check			
Гіте begin:		0704			Balance manufa	cturer	Mettler T	oledo
Гime end:		0705			Serial Number		C006982	938
Pass/Fail		Pass			Location of use		Chemistr	y lab
					Weight set seria	l# :	20410	
Pipet Calibrat	ion	0.01						
Dial setting: Time begin:	0655	0.01	units:	mL				
Replicate	Mass	Units		Volume (V)				
m <sub>1</sub>	0.0098	g		0.00983				
m <sub>2</sub>	0.0102	g		0.01024				
m <sub>3</sub>	0.0101			0.01014				
	0.0098	g		0.00983				
m <sub>4</sub>		g						
m <sub>5</sub>	0.0097	g		0.00973				
m <sub>6</sub>	0.0099	g		0.00993				
m <sub>7</sub>	0.0101	g		0.01014				
m <sub>8</sub>	0.0101	g		0.01014				
m <sub>9</sub>	0.0098	g		0.00983				
m <sub>10</sub>	0.0098	g		0.00983				
Mean mass	0.00993		PASS					
Гime end:	0705				COMMENTS:			
Mean vol.	0.0100				Mean mass accepta	ble range: ±8.0	0% or 0.00	08 g
3	0.0002				Mean mass acceptable range: 0.0092 to 0.0108 g			108 g
% CV	1.7794				CV tolerance accep	table range: <	3.0%	
CV Tolerance	PASS							
PASS/FAIL	PASS							
		6111						
					iter to the maximum nomi e, there should be no visib			
					mation is observed, succe			
	aken or the pip							
Complete green	-shaded cells o	nly (obtain	Z Corr. Factor from	Table 1 on Sheet 2). T	he spreadsheet will popul	ate all unshad	ed cells	
or you.								
Obtain RE and C	V tolerances fro	m the pipe	ttor manufacturer	's specifications.				
Jse appropriate	ASTM 1 or 2 tol	erances fo	r the balance chec	ks (refer to sheet 3).				



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 8 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

**Exhibit G11.3: Z-factor Correction Factor Table** 

	Air pressure, kPa (mbar)							
Temp degrees C	80 (800)	85 (850)	90 (900)	95 (950)	100 (1000)	101.3 (1013)	105 (1050)	
15	1.0017	1.0018	1.0019	1.0019	1.002	1.002	1.002	
15.5	1.0018	1.0019	1.0019	1.002	1.002	1.002	1.0021	
16	1.0019	1.002	1.002	1.0021	1.0021	1.0021	1.0022	
16.5	1.002	1.002	1.0021	1.0021	1.0022	1.0022	1.0022	
17	1.0021	1.0021	1.0022	1.0022	1.0023	1.0023	1.0023	
17.5	1.0022	1.0022	1.0023	1.0023	1.0024	1.0024	1.0024	
18	1.0022	1.0023	1.0023	1.0024	1.0025	1.0025	1.0025	
18.5	1.0023	1.0024	1.0024	1.0025	1.0025	1.0026	1.0026	
19	1.0024	1.0025	1.0025	1.0026	1.0026	1.0027	1.0027	
19.5	1.0025	1.0026	1.0026	1.0027	1.0027	1.0028	1.0028	
20	1.0026	1.0027	1.0027	1.0028	1.0028	1.0029	1.0029	
20.5	1.0027	1.0028	1.0028	1.0029	1.0029	1.003	1.003	
21	1.0028	1.0029	1.0029	1.003	1.0031	1.0031	1.0031	
21.5	1.003	1.003	1.0031	1.0031	1.0032	1.0032	1.0032	
22	1.0031	1.0031	1.0032	1.0032	1.0033	1.0033	1.0033	
22.5	1.0032	1.0032	1.0033	1.0033	1.0034	1.0034	1.0034	
23	1.0033	1.0033	1.0034	1.0034	1.0035	1.0035	1.0036	
23.5	1.0034	1.0035	1.0035	1.0036	1.0036	1.0036	1.0037	
24	1.0035	1.0036	1.0036	1.0037	1.0037	1.0038	1.0038	
24.5	1.0037	1.0037	1.0038	1.0038	1.0039	1.0039	1.0039	
25	1.0038	1.0038	1.0039	1.0039	1.004	1.004	1.004	
25.5	1.0039	1.004	1.004	1.0041	1.0041	1.0041	1.0042	
26	1.004	1.0041	1.0041	1.0042	1.0042	1.0043	1.0043	
26.5	1.0042	1.0042	1.0043	1.0043	1.0044	1.0044	1.0044	
27	1.0043	1.0044	1.0044	1.0045	1.0045	1.0045	1.0046	
27.5	1.0045	1.0045	1.0046	1.0046	1.0047	1.0047	1.0047	
28	1.0046	1.0046	1.0047	1.0047	1.0048	1.0048	1.0048	
28.5	1.0047	1.0048	1.0048	1.0049	1.0049	1.005	1.005	
29	1.0049	1.0049	1.005	1.005	1.0051	1.0051	1.0051	
29.5	1.005	1.0051	1.0051	1.0052	1.0052	1.0052	1.0053	
30	1.0052	1.0052	1.0053	1.0053	1.0054	1.0054	1.0054	

Z correction facotrs are for distilled water as a function of liquid test temperature and air pressure.

Reprinted from ISO 8655-6:2002(E).

If your laboratory does not have a mercury barometer, the true local barometric pressure for locations throughout North Carolina may be obtained from the North Carolina State University, State Climate Office of North Carolina website at http://www.nc-climate.ncsu.edu/.

Z values are in microliters per milligram.





SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 9 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

### Exhibit G11.4: Polypropylene Graduated Cylinder Volume Verification Logsheet.



Page 1 of 1

### Polypropylene Graduated Cylinder Volume Verification

# Pan balance calibration: Instrument Denver Instruments Company XL-3K Denver Instruments Company XL-3K Denver Instruments Company XL-3K Denver Instruments Company XL-3K Denver Instrument Services Certified Weights Certified weight (g) Certified weight (g) Certified weight (g)

Analyst initials	
	_

Analyst:		_	
Water and a	mbient te	moerature	

Cylinder type	Cylinder identification number	Volume Verified (mL)	Weight of volume in cylinder (g)	Corrected Volume* (mL) (weight of volume in cylinder x Z-factor)
100 mL	1	100 mL		
100 mL	2	100 mL		
100 mL	3	100 mL		
100 mL	4	100 mL		1
100 mL	5	100 mL		
100 mL	6	100 mL		
100 mL	7	100 mL		
100 mL	8	100 mL		
100 mL	9	100 mL		
100 mL	10	100 mL		
100 mL	11	100 mL		
100 mL	12	100 mL		
100 mL	13	100 mL		
100 mL	14	100 mL		
100 mL	15	100 mL		
100 mL	16	100 mL		
100 mL	17	100 mL		
100 mL	18	100 mL		
100 mL	19	100 mL		
100 mL	20	100 mL		
100 mL	21	100 mL		
100 mL	22	100 mL		
100 mL	23	100 mL		

<sup>\*</sup>The corrected volume for each cylinder must be within £5% of the true volume or volume verified (100 mL volume acceptance limit = 95 – 105 mL, if a cylinder exceeds these limits, it must be taken out of service.

SOP G11 - Revision 4 - Exhibit G11.4



SECTION SOP-G11
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 10 OF 10

## Subject: Mechanical Pipette Calibration and Polypropylene Cylinder Volume Verification

Exhibit G11.5: Polypropylene Graduated Cylinder Volume Verification excel spreadsheet

	Polyp	ropylene Grad	luated Cylin	der Volume V	erification
Pan balance calibrat	ion:				
Instrument	Serial number	Calibration company			
Denver Instruments Company XL-3K	B039128	Laboratory Instrument Services			
Certified Weights	20410	Precision Weighting			
Date:	12-17-24				
Analyst:	JP				
Water and ambient temperature:	23.3°C				
Z Factor:	1.0036				
	Cylinder type	Cylinder identification number	Volume Verified (mL)	Weight of volume in cylinder (g)	Corrected Volume* (mL) (weight of volume in cylinder x Z-factor)
	100 mL	1	100 mL	98	98
	100 mL	2	100 mL	97	97
	100 mL	3	100 mL	97	97
	100 mL	4	100 mL	99	99
	100 mL	5	100 mL		
	100 mL	6	100 mL	99	99
	100 mL	7	100 mL	98	98
	100 mL	8	100 mL	98	98
	100 mL	9	100 mL	99	99
	100 mL	10	100 mL	99	99
	100 mL	11	100 mL	98	98
	100 mL	12	100 mL	98	98
	100 mL	13	100 mL	98	98
	100 mL	14	100 mL	99	99
	100 mL	15	100 mL	97	97
	100 mL	16	100 mL	98	98
	100 mL	17	100 mL	99	99
	100 mL	18	100 mL	97	97
	100 mL	19	100 mL	99	99
	100 mL	20	100 mL	99	99
	100 mL	21	100 mL	99	99
	100 mL	22	100 mL	100	100
	100 mL	23	100 mL	97	97
	100 mL	24	100 mL	97	97



SECTION SOP-G12
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 1 OF 5

## **Subject: Thermometer Calibration**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Type		
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
09-01-09	1	External	William Rogers (TVA)	Incubator, refrigerator, and drying oven temperatures section and
		(TVA,	Cynthia Russell (TVA)	associated exhibits removed and added to SOP C1.
		Environmental	Rick Sherrard (TVA)	
		Standard, Inc.)	Rock Vitale	
			(Environmental	
			Standards, Inc.)	
		Internal	Jim Sumner (ETS)	
06-01-25	2	Internal	Jim Sumner (ETS)	Updated reference to 24 <sup>th</sup> Edition of Standard Methods.

## **Scope and Application**

To calibrate thermometers used to measure the temperature of water samples used in toxicity tests, wastewater, receiving water, drinking water, and incubators/refrigerators.

## **Quality Control**

**Standardization**: All mercury, red spirit-filled, hand-held thermometers, and meters, which measure temperature, must be verified at least <u>annually</u> (once every calendar year) with traceable NIST thermometers.

In addition, any thermometer used to monitor the temperature of incubators or water baths must be verified at least **quarterly** with traceable NIST thermometers.

All digital thermometers must be verified at least **quarterly** with traceable NIST thermometers.



SECTION SOP-G12
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 2 OF 5

## **Subject: Thermometer Calibration**

### **Equipment and Materials**

Mercury-filled thermometers
Hand-held thermometers
Digital thermometers
NIST traceable thermometers
Rinse bottle
Deionized water
Waste container
Sand
Thermometer Calibration Logsheet

### **Procedure**

#### A. Calibration.

- All mercury-filled and hand-held thermometers and meters, which measure temperature, must be calibrated with a traceable NIST thermometer (refer to Exhibit G12.1 for NIST certification records).
  - a. Thermometers used to measure the temperature within temperature-controlled incubators are calibrated at  $20.0 \pm 1.0^{\circ}$ C,  $25.0 \pm 1.0^{\circ}$ C, or  $35.0 \pm 0.5^{\circ}$ C.
  - b. Thermometers used to measure the temperature within refrigerators are calibrated at 0.1  $6.0^{\circ}$ C.
  - c. Thermometers used to measure the temperature within drying ovens are calibrated at  $60.0 \pm 2.0$  °C,  $105 \pm 2.0$  °C, or  $180 \pm 2.0$  °C.
  - d. General use thermometers are calibrated at 25.0 ± 1.0 °C and 0.1 4.0 °C.
- 2. Place all thermometers to be calibrated at the desired temperature in a beaker of water (or sand for ovens) with the NIST thermometer. Place the beaker of water (or sand) in the respective temperature-controlled incubator for calibrating thermometers at  $20.0 \pm 1.0^{\circ}$ C,  $25.0 \pm 1.0^{\circ}$ C, or  $35.0 \pm 0.5^{\circ}$ C, a refrigerator for calibrating thermometers at  $0.1 \pm 6.0^{\circ}$ C, or the respective drying oven for calibrating thermometers at  $60.0 \pm 2.0^{\circ}$ C,  $105 \pm 2.0^{\circ}$ C, or  $180 \pm 2.0^{\circ}$ C. Allow the water (or sand) to reach the desired temperature. For each thermometer, compare the temperature measurement to the measurement on the NIST thermometer. If the thermometer has an immersion line, it must be utilized. If



SECTION SOP-G12
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 3 OF 5

## **Subject: Thermometer Calibration**

it does not have an immersion line, both calibration and measurements must be made under the same water depth.

- 3. Record the date of calibration, serial number of the thermometer, calibration temperature, and correction factor (temperature adjustment of the thermometer to equal the measurement of the NIST thermometer) on the Thermometer Calibration Logsheet (Exhibit G12.2) and on a label affixed to the thermometer.
- 4. Thermometers having breaks in the mercury or having a correction factor > 1.0 are taken out of service.

### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

### **Exhibits**

Exhibit G12.1: NIST Thermometer Certification Records. Exhibit G12.2: Example Thermometer Calibration Logsheet.



SECTION SOP-G12
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 4 OF 5

## **Subject: Thermometer Calibration**

#### **Exhibit G12.1: NIST Thermometer Certification Records.**



### **Precision Weighing**

1949 Evans Road Cary, North Carolina 27513 Phone: (919) 678-0077 \* Fax: (919) 678-0078 Email: pweighing@aol.com

Client: Contact: Department: ETS Inc. Jim Sumner Lab

Calibration SOP: Description: Digital Thermometer 2020 Rev 1.1 Instrument ID: 61786906 Serial Number: 61786906 Manufacturer: Fisher Scientific Operating Range: 0 to 180 °C Model: 15-078-2 Instrument Range: -350 to 2200 °C Calibration Interval: Annually Calibration Tolerance: +/- 1 deg C

Calibration Notes: Nominal Targets are 0,20,25,35,44.5,60,105,150,180

Calibration Data ("As left" data is identical to "As found" data, if "As left" is blank)

Test Points	Units	Standard	As Found	Error	As Left	Error	Pass/Fail
0	°C	0.0036	0.1	0.0964			Pass
20	° C	20.0022	19.9	-0.1022			Pass
25	° C	24.9429	24.9	-0.0429			Pass
35	°C	34.9700	34.8	-0.1700			Pass
44.5	* C	44.3375	44.4	0.0625			Pass
60	° C	59.8021	59.8	-0.0021			Pass
105	° C	104.690	105.0	0.310			Pass
150	° C	149.333	150.0	0.667			Pass
180	° C	179.713	180.2	0.487			Pass

Instrument Found in Tolerance? Yes
Instrument Left in Tolerance? Yes

Calibration Date:	20 Nov 24	Calibration Due Date:	Nov 25
Test Standards:	Standard ID:	Expiration Date:	
Digital Thermometer	230670827	03 Oct 25	
Technicians Remarks:	None.		
Calibration Performed By		Approved Signature:	
	20 Nov 24	Jim Sumne	





SECTION SOP-G12
REVISION NUMBER 2
EFFECTIVE DATE 06-01-25
PAGE 5 OF 5

## **Subject: Thermometer Calibration**

### **Exhibit G12.2: Example Thermometer Calibration Logsheet.**



Page 1 of 5

#### **Thermometer Calibration**

Traceable Thermometer:	Fisher Scientific Digital Thermometer	Date:	
Probe Serial number:	61786906	Analyst:	
Meter Serial number:	61789314		
Calibration date:	Probe and Meter: 11-20-24	Ice Point of Traceable Thermometer:	
Re-calibration due:	Probe and Meter: 11-2025		

Thermometer Serial Number	Location	Location Serial Number	Thermometer Temperature (°C)	Traceable Thermometer Temperature (°C)	Correction Factor (°C)
	Refrigerator # 1	BA64722418	Not in use.		
96-01587	Refrigerator # 2	BA65013194			
96-01566	Refrigerator # 3	BA11411547			
6959	Refrigerator # 4 TOP	WA42501619			
4934	Refrigerator # 4 BOTTOM	WA42501619			
4565	Refrigerator #5 TOP	WA42500397			
6315	Refrigerator # 5 BOTTOM	WA42500397			
95-02120	Refrigerator – Algae/YWT	E2001648367			
WB72043714	BOD Incubator # 1	WB72043714	Calibrated to 20.0°C		
7126	BOD Incubator # 1 Thermometer	WB72043714			
WB94157442	BOD Incubator # 2	WB94157442	Calibrated to 20.0°C		
7049	BOD Incubator # 2 Thermometer	WB94157442			
WB22114719	BOD Incubator # 3.	WB94157442	Calibrated to 20.0°C		
6209	BOD Incubator # 3 Thermometer	WB94157442			
WB42676093	Reagent incubator # 1	WB42676093	Calibrated to 25.0°C		
5030	Reagent Incubator # 1	WB42676093			
WB41340472	Toxicity Incubator # 1	WB41340472	Calibrated to 25.0℃		
160724968	Toxicity Incubator # 1 MIN/MAX	WB41340472			
6272	Toxicity Incubator # 1 TOP	WB41340472			
4673	Toxicity Incubator # 1 BOTTOM	WB41340472			
WB22114693	Toxicity Incubator # 2	WB22114693	Calibrated to 25.0℃		
130761061	Toxicity Incubator # 2 MIN/MAX	WB22114693			
4676	Toxicity Incubator # 2 TOP	WB22114693			
4755	Toxicity Incubator # 2 BOTTOM	WB22114693			
WB95219633	Toxicity Incubator # 4	WB95219633	Calibrated to 25.0°C		
160761060	Toxicity Incubator # 4 MIN/MAX	WB95219633			
4547	Toxicity Incubator # 4 TOP	WB95219633			
6216	Toxicity Incubator # 4 BOTTOM	WB95219633			

Note: If correction factor is ≥ 1.0°C, the thermometer must be taken out of service

SOP G12-Revision 2-Exhibit G12.2



SECTION SOP-G13
REVISION NUMBER 2
DATE 06-01-25
PAGE 1 OF 3

## **Subject: Desiccator Maintenance**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Julunoe-	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Туре		
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
03-01-20	1	External (TVA)	Rick Sherrard (TVA)	Updated procedure and benchsheet to include the serial number of the meter used to perform chlorine.
		Internal	Jim Sumner (ETS)	
06-01-25	2	Internal	Jim Sumner (ETS)	Updated procedure during document review.

## **Scope and Application**

To ensure proper desiccator maintenance by drying the desiccant.

## **Quality Control**

Desiccator is to be dried at least monthly.

## **Equipment and Materials**

Desiccator with plastic desiccant pans Color Indicating Desiccant Desiccator Maintenance Log 103 -105 °C forced air oven Metal drying pan High temperature gloves



SECTION SOP-G13
REVISION NUMBER 2
DATE 06-01-25
PAGE 2 OF 3

## **Subject: Desiccator Maintenance**

### **Procedure**

### A. Drying the Desiccant

- 1. Empty the desiccant into a metal drying pan. Desiccant should be blue if dry or purple if it has absorbed moisture.
- 2. Place the metal pan into the 103 105°C oven until desiccant turns blue.
- 3. After the desiccant turns blue, use the high temperature gloves and pour the desiccant into the plastic desiccant trays.
- 4. Place the plastic trays into the bottom of each desiccator.
- 5. Record the date and the desiccant is dried in the desiccator maintenance log.
- 6. Desiccator is to be dried at least monthly.

#### **Exhibits**

Exhibit G13.1: Desiccator Maintenance Benchsheet.



SECTION SOP-G13
REVISION NUMBER 2
DATE 06-01-25
PAGE 3 OF 3

## **Subject: Desiccator Maintenance**

### Exhibit G13.1: Desiccator Maintenance Benchsheet.



### Page \_\_\_\_

### **Desiccant Logsheet**

Desiccator #1	Dessicator #2	Dessicator #3	Dessicator #4	Analyst
	Desiccator #1	Desiccator #1 Dessicator #2	Desiccator #1 Dessicator #2 Dessicator #3	Desiccator #1 Dessicator #2 Dessicator #3 Dessicator #4

SOP G13-Revision 2 - Exhibit G13.1



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 1 OF 6

**Subject: Illumination** 

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

### **Document Revision History**

Effective	Revision	Review Type	Evaluators	Revisions
Date	number			
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
01-03-12	1	Internal	Jim Sumner (ETS)	Updated exhibits during document review.
03-01-20	2	External (TVA)	Rick Sherrard (TVA)	Updated procedure and bench sheet to describe the current quality control and measurement of luminosity.
		Internal	Jim Sumner (ETS)	
06-01-25	3	Internal	Jim Sumner (ETS)	Updated exhibits during document review.

## **Scope and Application**

To document the photoperiod and intensity of light necessary for maintaining cultures and toxicity tests.

## **Quality Control**

- 1. The intensity of light (luminosity) in all incubators used for maintaining cultures and toxicity tests and the algae culture area must be measured at least **quarterly** or when light fixtures are replaced. The luminosity should be maintained between 50 100 foot-candles (ft-c) for maintaining cultures and toxicity tests and between 360 440 foot-candles for maintaining algae cultures. For toxicity tests in support of AL NPDES discharges, the luminosity of incubators should be measured each day that the tests are performed. This is documented on the test specific bench sheets.
- 2. A uniform photoperiod of 16-hours light and 8-hours dark must be maintained in all incubators used for maintaining cultures and toxicity tests as well in the fish culture laboratory. Automatic timers and internal timers within the incubators regulate this photoperiod. The correct time displayed on these timers is documented and adjusted according to SOP-G16. The correct



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 2 OF 6

### **Subject: Illumination**

photoperiod is documented **quarterly** when luminosity is documented. Timers may need to be adjusted when incubators are serviced. Algae cultures are maintained under constant illumination.

## **Equipment and Materials**

Light meters
Incubators
Light fixtures
Window screen or black electrical tape
Automatic timers
Luminosity Log Sheet

Procedure (Meter: TACKlife LM01 Light Meter, SN 201612033203)

### A. Verification of Light Meter.

Prior to taking luminosity measurements, verify the luminosity of the meter. To verify the Tacklife LM01 Light Meter, compare the luminosity displayed to a separate meter (Leaton L860 LUX Meter, SN not provided). Place each probe face up on the surface of an area to be measured. Each display must be with 10% of one another prior to taking luminosity measurements. Measure and record the luminosity of each meter and percent difference on the Luminosity Log Sheet (Exhibit G14.1).

### B. Measuring Luminosity.

- 1. Turn the light meter on and set the meter to read in foot-candles (ft-c).
- Place the probe with the sensor face up on the surface of the area to be measured. Check the luminosity throughout the area where cultures or tests are maintained. In the fish culture laboratory to monitor the luminosity of each tank, hold the probe with the sensor face up close to the water surface.
- 3. If the luminosity is outside the 50 100 ft-c suggested range (360 440 ft-c for algae), adjust the light intensity by adding or removing light fixtures. Window screen or black electrical tape may also be placed over light fixtures to adjust the light intensity in specific areas.
- 4. Once the luminosity has been corrected in all areas, measure and record the luminosity at the four corners and middle of each shelf in each incubator, the location were algae is



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 3 OF 6

### **Subject: Illumination**

cultured and at the water surface for each tank in the fish culture laboratory. Allow readings to stabilize before recording in the Luminosity Log. Measurements must be in foot-candles (ft-c) and reported to the nearest 1 ft-c.

5. Verify and document in the Luminosity Log that the time is correct, and the photoperiod program is set for 16-hour light and 8-hours dark for each incubator and in the fish culture laboratory.

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

### References

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms, 4<sup>th</sup> ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-014. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.

**Instrument Manuals** 

#### **Exhibits**

Exhibit G14.1: Luminosity Log Sheet.



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 4 OF 6

## **Subject: Illumination**

### **Exhibit G14.1: Luminosity Log Sheet.**



Page 1 of 2

### **Luminosity Log**

Meter:	TACKlife Digital LUX Meter, LM01 201612033203		Date:	
Serial number:			Analyst:	
Purchase date:	06-23-17	06-23-17		
Light Meter	Comparison to Leaton L830	TACKlife (ft-c)		
Verification	Lux Meter SN not provided. (must be ± 10% prior to taking measurements)	Leaton (ft-c)		
		% Difference		

Incubator number	Shelf	Luminosity (ft-c) at location indicated on shelf Acceptable range = 50 – 100 ft-c	Incubator number	Shelf	Luminosity (ft-c) at location indicated on shelf Acceptable range = 50 – 100 ft-c
Incubator 1 SN: WB41340472	Α	< 50 Do not use.	Incubator 2 SN: WB22114693	Α	< 50 Do not use.
	В			В	
	с			с	
	D			D	
	E			E	
	F			F	
	Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.		Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.
Incubator 4 SN: WB95219633	Α	< 50 Do not use.	Incubator 5 SN: WB60643874	Α	< 50 Do not use.
	В			В	
	С			с	
	D			D	
	E			E	
	F			F	
	Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.		Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.

SOP G14-Revision 3- Exhibit G14.1



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 5 OF 6

## **Subject: Illumination**



Page 2 of 2

### **Luminosity Log**

Meter:	TACKlife Digital LUX Meter, LM01	Date:	- 1
Serial number:	201612033203	Analyst:	

Incubator number	Shelf	Luminosity (ft-c) at location indicated on shelf Acceptable range = 50 – 100 ft-c	Incubator number	Shelf	Luminosity (ft-c) at location indicated on shelf Acceptable range = 50 – 100 ft-c.	
Incubator 6 SN: WB42667925	Α	< 50 Do not use.	Incubator 7 SN: WB42668003	Α	< 50 Do not use.	
	В			В		
	c			c		
	D			D		
	E			E		
	F			F		
	Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.		Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.	
	Α	< 50 Do not use.	Algae Culture Area	Luminosity (ft-c) at location indicated Acceptable range = 360 – 440 ft-c		
	В					
	С					
	D					
	E					
	F					
	Initials	Verified photoperiod program (16-hours light and 8-hours dark) and current time.				

SOP G14-Revision 3- Exhibit G14.1



SECTION SOP-G14
REVISION NUMBER 3
EFFECTIVE DATE 06-01-25
PAGE 6 OF 6

**Subject: Illumination** 



Page 1 of 1

#### Luminosity Log, Fish Culture Laboratory

Meter:	TACKlife Digital LUX Me	Date:		
Serial number:	201612033203		Analyst:	
Purchase date:	06-23-17			
Light Meter	Comparison to Leaton L830 Lux	TACKlife (ft-c)		
Verification	Meter 5N not provided. (must be ± 10% prior to taking	Leaton (ft-c)		
	measurements)	% Difference		

Verified photoperiod program (16-hours light and 8-hours dark)	Analyst:	
and current time.	150-000	

Breeding Tank Set	ID#	Luminosity (ft-c) at water surface Acceptable range = 50 – 100 ft-c	Breeding Tank Set	ID#	Luminosity (ft-c) at water surface Acceptable range = 50 – 100 ft-c
Α	1		В	1	
	2			2	
	3		7 0	3	
	4		1 1	4	
	5			5	
	6			6	
С	1		D	1	
- X-	2			2	
	3			3	
	4			4	
	5			5	
	6			6	
E	1		F	1	
-	2			2	
	3			3	
	4			4	
	5		311 A	5	
	6		3	6	
G	1		Н	1	
	2			2	
	3		1 1	3	
	4		7 - 3	4	
	5			5	
	6			6	

Stock tank 1: \_\_\_\_ 2: \_\_\_ 3: \_\_\_ 5: \_\_\_ 6: \_\_\_ 7: \_\_\_ 8: \_\_\_ SOP G14-Revision 3- Exhibit G14.1



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 1 OF 10

## **Subject: Reagent, Stock Standard, and Chemical Preparation**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	~	06-01-25
Quality Assurance Officer	Jim Sumner	Jun/unse-	06-01-25

## **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Type		
12-01-00	0	Internal	Jim Sumner (ETS)	Original document
09-01-09	1	External	William Rogers (TVA)	Updated exhibits during document review.
		(TVA,	Cynthia Russell (TVA)	Exhibits modified to include the initials and date of the author making
		Environmental	Rick Sherrard (TVA)	entries to the logbooks. Logbooks modified to provide additional
		Standard, Inc.)	Rock Vitale	traceability.
			(Environmental	Provided further clarification on the use of expiration dates.
			Standards, Inc.)	
		Internal	Jim Sumner (ETS)	
01-03-12	2	Internal	Jim Sumner (ETS)	Updated exhibits during document review.
				Provided further clarification on the use of expiration dates.
06-20-12	3	External	William Rogers (TVA)	• Expiration dates changed from 1 year or 5 years from open to 1 year or 5
		(TVA)	Donald Snodgrass (TVA)	years from <u>receipt</u> .
			Rick Sherrard (TVA)	
		Internal	Jim Sumner (ETS)	
03-01-20	4	External	Rick Sherrard (TVA)	Provided additional guidance for chemicals and consumables used for
		(TVA)		toxicity testing, which require testing prior to use.
		Internal	Jim Sumner (ETS)	
06-01-25	5	Internal	Jim Sumner (ETS)	Updated reference to 24 <sup>th</sup> Edition of Standard Methods.

# **Scope and Application**

To provide traceability of all reagents, stock standards, and chemicals used in the laboratory.

# **Toxicity Testing Quality Control**

In addition to being "checked in" (as described in this SOP), analytical testing or toxicity "checks" may be required of certain chemicals or consumables used for toxicity testing. These chemicals or consumables must be tagged as "do not use" prior to testing. After testing is performed and the chemicals or consumables meet quality control requirements, these items will be tagged with "Acceptable to Use" and the date that testing was performed. At this time, these chemicals or consumables are cleared for



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 2 OF 10

## **Subject: Reagent, Stock Standard, and Chemical Preparation**

use. Chemicals and consumables that fall under this requirement include: cubitainers, medicine cups, Solo cups, MarineMix and synthetic water chemicals used in toxicity tests.

### **Equipment and Materials**

Reagent log Stock standard log Chemical log

#### **Procedure**

- A. "Checking-In" Reagents, Stock Standards, Chemicals and Consumables.
  - 1. All reagents, stock standards, and chemicals must be checked-in and assigned an identification number, date received, date opened sticker (or placed into use), and a safety label. Information is recorded in one of the following logs (Exhibit G15.1):
    - Chemical Log
    - Reagent Log
    - Stock Standard Log
  - 2. Chemicals, stock standards, and reagents are used on or before their expiration dates. At the time of analysis, any unused, expired chemicals, stock standards, and reagents are disposed. Expiration dates for chemicals, stock standards, and reagents (that are not assigned as expiration date by the manufacturer) are assigned in a MONTH DAY YEAR format. These chemicals, stock standards, and reagents expire on the expiration date.
  - 3. The chemical log is used to record information for all chemicals and consumables (e.g. Petri dishes, filters). Each chemical or consumable is assigned an identification number (CHM#).
    - a. Each chemical receives a label with the CHM#, date received, expiration date, date opened (or placed into use), and initials. The expiration date for dry chemicals that are not assigned an expiration date by the manufacturer is five years from the date received. The expiration date for wet chemicals, that are not assigned an expiration date by the manufacturer, is one year from the date received. Consumables, such as containers, medicine cups, solo cups, acetone or acid used for cleaning purposes only, are not assigned an expiration date.



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 3 OF 10

- b. If a Safety Data Sheet (SDS) is not received with the chemical or is not contained in the SDS Log, contact the supplier to obtain the SDS. SDS's for each chemical must be stored in the SDS Log. As new chemicals are received, SDS's are replaced with the most current SDS received from the manufacturer.
- c. If the chemical is not compatible with other types of chemicals, it must be stored in a separate location. This may include acids, oxidizers, ammonium hydroxide, flammables, etc.
- d. If more than one container of the same chemical is received, assign each container the same chemical number (provided they are the same lot or batch). For example. If two containers are received, they should be labeled "1 of 2" and "2 of 2".
- e. Record the following information for the chemical or consumable in the chemical log:
  - Page number
  - Chemical number
  - Analyst
  - Date
  - Chemical name
  - Test or application used for
  - Material type
  - Toxicity Testing Quality Control required?
  - Manufacturer
  - Lot number
  - Chemical concentration
  - Number or volume received
  - Date received
  - Expiration date
  - Received by
  - Comments
  - Storage
  - Hazards
  - Personal hygiene (health, fire, reactivity, personal equipment)



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 4 OF 10

- 4. The reagent log is used to record information for all reagents. Each reagent is assigned an identification number (INR#).
  - a. Each reagent receives a label with the INR#, date prepared or received, expiration date (1 year from preparation or manufacturer expiration date), date opened (or placed into use) and initials. The expiration date for reagents, that are not assigned an expiration date by the manufacturer, is one year from the date received or preparation date.
  - b. If an SDS is not received with the reagent or is not contained in the SDS Log, contact the supplier to obtain the SDS. SDS's for each reagent must be stored in the SDS Log. As new reagents are received, SDS's are replaced with the most current SDS received from the manufacturer.
  - c. If the reagent is not compatible with other types of chemicals, it must be stored in a separate location. This may include acids, oxidizers, ammonium hydroxide, flammables, etc.
  - d. If more than one container of the same reagent is prepared or received, assign each container the same reagent number (provided they are the same lot or batch). For example. If two containers are prepared or received, they should be labeled "1 of 2" and "2 of 2".



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 5 OF 10

- e. Record the following information for the reagent in the reagent log:
  - Page number
  - Reagent number
  - Analyst
  - Date
  - Reagent name
  - Test or application used for
  - Type
  - Manufacturer
  - Lot number
  - Reagent concentration
  - Number and/or volume
  - Date received/prepared
  - Expiration date
  - Received/prepared by
  - Directions for preparing the reagent
  - Reagent prepared from
  - Storage
  - Hazards
  - Personal hygiene (health, fire, reactivity, personal equipment)
- 5. The stock standard log is used to record information for all stock standards. Each stock standard is assigned an identification number (INSS#).
  - a. Each stock standard receives a label with the INSS#, date prepared or received, expiration date (1 year from preparation/received or manufacturer expiration date), date opened (or placed into use) and initials. The expiration date for stock standards, that are not assigned an expiration date by the manufacturer, is one year from the date received or preparation date.
  - b. If the stock standard is not compatible with other types of chemicals, it must be stored in a separate location. This may include acids, oxidizers, ammonium hydroxide, flammables, etc.
  - c. If more than one container of the same stock standard is prepared or received, assign each container the same stock standard number (provided they are the same lot or batch). For example. If two containers are prepared or received, they should be labeled "1 of 2" and "2 of 2".



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 6 OF 10

- f. Record the following information for the stock standard in the stock standard log:
  - Page number
  - Stock standard number
  - Analyst
  - Date
  - Standard name
  - Test or application used for
  - Type
  - Manufacturer
  - Lot number
  - Stock concentration
  - Number and/or volume
  - Date received/prepared
  - Expiration date
  - Received/prepared by
  - Directions for preparing the stock
  - Stock prepared from
  - Storage
  - Hazards
  - Personal hygiene (health, fire, reactivity, personal equipment)



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 7 OF 10

## **Subject: Reagent, Stock Standard, and Chemical Preparation**

### **Safety and Hazardous Waste Management**

Safety glasses, gloves and lab coats should always be worn.

Review Policy-P6: General Safety Policy and Policy-P9: Radiation Protection Policy for additional safety requirements.

#### References

Standard Methods for the Examination of Water and Wastewater, 24<sup>th</sup> Edition, 2023. American Public Health Association, 800 I Street, NW, Washington DC 20001-3710.

USEPA. October 2002. Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms, 5<sup>th</sup> ed. EPA-821-R-02-012. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms, 4<sup>th</sup> ed. EPA-821-R-02-013. US Environmental Protection Agency, Cincinnati, OH.

USEPA. October 2002. Short-Term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Marine and Estuarine Organisms, 3<sup>rd</sup> ed. EPA-821-R-02-014. US Environmental Protection Agency, Cincinnati, OH.

TNI Standard. Management and Technical Requirements for Laboratories Performing Environmental Analysis. EL-V1-ISO-2016-Rev2.0. The NELAC Institute, PO Box 2439, Weatherford, TX 76086.

#### **Exhibits**

G15.1: Example Chemical, Reagent and Stock Standard Logs



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 8 OF 10

# **Subject: Reagent, Stock Standard, and Chemical Preparation**

## **G15.1:** Example Chemical, Reagent and Stock Standard Logs

Invironmental Testing Solutions, Inc.	The state of the s
and the same teat.	Chemical Log
hemical number: CHM_	
Analyst:	
Date:	
Chemical name:	
Test or application used for:	
Material type:	□ Consumable □ Chemical
Toxicity Testing:	
Quality Control required?	□ Yes □ No
Manufacturer:	<u> </u>
Lot number:	
Chemical concentration:	
Number and/or volume:	
Date received:	
Expiration date:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on	Wet Chemicals = 1 year from receipt, if not assigned by manufacture
Expiration date:  The Chemical identified on  Received by:	Dry Chemicals = 5 years from receipt, if not assigned by manufacture Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture this form was discarded on or before the expiration date listed above.
Expiration date:  The Chemical identified on  Received by:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on  Received by:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture  Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on  Received by:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture  Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on  Received by:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on  Received by:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on  Received by:  Comments:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:  Storage: Hazards:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture
Expiration date:  The Chemical identified on Received by: Comments:  Storage: Hazards:	Wet Chemicals = 1 year from receipt, if not assigned by manufacture Consumables = no expiration date, unless assigned by manufacture this form was discarded on or before the expiration date listed above.

30F G13-NEVISION 3-EXHIBIT G13.1



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 9 OF 10

# **Subject: Reagent, Stock Standard, and Chemical Preparation**

CTC				
Environmental Testing Solution	is, Inc.	Et. 80.		
	4.5	Reagent I	.og	
eagent num	ber: INR			
Analyst:				
Date:				
leagent name				
est or applica	tion used for:			
ype:		☐ Purchased ☐ Laborato	ry prepared	
Manufacturer:				
ot number:				
Reagent conce	ntration:			
Number and/o	r volume:			
Number and/or volume:				
Date received/	prepared:			
Expiration date	e:	(1 ye		if not assigned by manufacturer date listed above.
Date received/ Expiration date The Reage Received/prep	e: ent identified on t			
The Reage	e: ent identified on t	nis form was discarded on or b		
The Reage Received/prep Directions for	e: int identified on the lared by: preparing reago	nis form was discarded on or b		
The Reage Received/prep Directions for p	e: int identified on the lared by: preparing reago	nis form was discarded on or b		date listed above.
The Reage Received/prep Directions for p	e: int identified on the lared by: preparing reago	nis form was discarded on or b		date listed above.
The Reage Received/prep Directions for p	e: int identified on the lared by: preparing reago	nis form was discarded on or b		date listed above.
The Reage Received/prep Directions for p Reagent prepa Storage:	e: int identified on the lared by: preparing reago	nis form was discarded on or b		date listed above.
The Reage	e: int identified on the sared by: preparing reagons red from:	nis form was discarded on or b		date listed above.

SOP G15-Revision 5- Exhibit G15.1



SECTION SOP-G15
REVISION NUMBER 5
EFFECTIVE DATE 06-01-25
PAGE 10 OF 10

# **Subject: Reagent, Stock Standard, and Chemical Preparation**

	3			
Environmental Testing 1	Solutions, Inc.			
		Stock Standard Log		
ock stand	dard number: I			
out stant	adra Hambert II			
Analyst:				
Date:				
Sec. 5 5 -				
tandard na				
est or app	lication used for:			
Гуре:		☐ Purchased ☐ Laboratory prepared		
Manufactu	rer:			
ot number	•			
stock conce	entration:			
Number an	d/or volume:			
Date receiv	ed/prepared:			
Expiration date:				
		(1 year from preparation/re on this form was discarded on or before the exp	ceipt, if not assigned by manufacturer iration date listed above.	
The Stock				
The Stock	Standard identified	on this form was discarded on or before the exp		
The Stock	Standard identified	on this form was discarded on or before the exp		
The Stock	Standard identified	on this form was discarded on or before the exp		
The Stock	Standard identified	on this form was discarded on or before the exp		
The Stock Received/p Directions f	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f	Standard identified	on this form was discarded on or before the exp		
The Stock Received/p Directions f	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f Stock stand Storage:	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f Stock stand Storage: Hazards:	Standard identified repared by: for preparing stock	on this form was discarded on or before the exp	iration date listed above.	
The Stock Received/p Directions f	Standard identified repared by: for preparing stock lard prepared from Health:	this form was discarded on or before the exp	iration date listed above.	

SOP G15-Revision 5- Exhibit G15.1



SECTION SOP-G16
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 1 OF 3

## **Subject: Time Calibration**

## **Approval**

Title	Name	Signature	Date
Laboratory Supervisor	Kelley E. Keenan	N	06-01-25
Quality Assurance Officer	Jim Sumner	Infune	06-01-25

# **Document Revision History**

Effective	Revision	Review	Evaluators	Revisions
Date	number	Type		
06-01-11	0	Internal	Jim Sumner (ETS)	Original document
06-01-11	1	External	Cynthia Russell (TVA)	Method for calibrating, recording and documenting time.
		(TVA)	Rick Sherrard (TVA)	
		Internal	Jim Sumner (ETS)	
09-01-11	2	External	Rick Sherrard (TVA)	Corrected typographical error.
		(TVA)		
		Internal	Jim Sumner (ETS)	
03-01-20	3	Internal	Jim Sumner (ETS)	Updated procedure to include NELAP requirements.
06-01-25	4	Internal	Jim Sumner (ETS)	Updated procedure to calibrate clocks to NIST Official U.S. Time.

# **Scope and Application**

To calibrate various clocks used to record and document time throughout the laboratory

# **Quality Control**

All clocks, which are used to document or record time, must be calibrated quarterly.

# **Equipment and Materials**

Various clocks NIST Official U.S. Time Cell phone Time Calibration Log



SECTION SOP-G16
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 2 OF 3

## **Subject: Time Calibration**

#### **Procedure**

- 1. All clocks, which are used to document or record time, must be calibrated quarterly.
- 2. Clocks are calibrated to the nearest minute, since the recording of time on individual log sheets is documented to the nearest minute in the laboratory.
- 3. The NIST Official U.S. Time is used as the standard to calibrate clocks within the laboratory. https://time.gov
- 4. A cell phone is calibrated to the NIST Official U.S. Time as indicated above and will be used to calibrate all individual clocks throughout the laboratory. This is documented in Time Calibration Log (Exhibit G16.1).

#### References

Time and Frequency Division, National Institute of Standards and Technology (NIST) Measurement Laboratory, US Department of Commerce.

Instrument Manual

#### **Exhibits**

Exhibit G16.1: Time Calibration Log.



SECTION SOP-G16
REVISION NUMBER 4
EFFECTIVE DATE 06-01-25
PAGE 3 OF 3

**Subject: Time Calibration** 

#### Exhibit G16.1: Time Calibration Log.



Page\_\_\_\_

#### **Time Calibration Log**

Standard:

NIST Official U.S. Time

https://time.gov

Date	NIST U.S. Official Time	All Clocks Calibrated to NIST Official U.S. Time (Analyst Initials)
	1	

SOP G16-Revision 4-Exhibit G16.1